



Corporate Headquarters
400 Valley Road
Warrington, PA 18976
1-800-523-2575
FAX 1-800-343-3291
Email: info@polysciences.com
www.polysciences.com

Europe - Germany
Polysciences Europe GmbH
Handelsstr. 3
D-69214 Eppelheim, Germany
(49) 6221-765767
FAX (49) 6221-764620
Email: info@polysciences.de

TECHNICAL DATA SHEET 289

Page 1 of 1

Tissue-Tack

Adhesive for JB-4[®], JB-4[®] Plus and Immunol/Bed

Background:

Loss of sections during lengthy histochemical procedures causes frustration and diagnostic delays. Tissue-Tack is an acrylic-based adhesive which has been used successfully to attach JB-4[®] sections to slides when conventional section adhesives have proved to be insufficient.

Procedure:

Application of Tissue-Tack in the following method allows for the effective demonstration of immunofluorescent and immunoperoxidase staining. This adhesive is also applicable to the Esterbauer and Nöhammer modification of the DDD test.

1. Dry microtomed JB-4[®] sections onto a clean glass slide.
2. A small perforation should be made in the JB-4[®] plastic around the perimeter of the tissue section, to expose the glass. Searing the edges of the plastic with a hot instrument provides additional mechanical adhesion.
3. Using an applicator stick, apply a small amount of Tissue-Tack, covering the perforated area of the JB-4[®] section.
4. Quickly smear the adhesive in the direction away from the section to form a thin film over the slide. (This procedure should be performed under a fume hood.)
5. Air dry for a few minutes in the fume hood, until the adhesive is clear and unsticky.
6. Proceed with desired procedures. Tissue-Tack becomes soft and tacky when it is subjected to strong alkalis and solvents, and adheres to the glass slide, but not to the JB-4[®]. It can be removed after the staining techniques.

Caution:

Prolonged exposure can cause irritation to eyes, skin and respiratory tract. Wear protective glasses and gloves. Use only with adequate ventilation. Keep container tightly closed.

Reference:

Pappas, Charles N., Stain Tech., 57, 357 (1982).

Ordering Information:

Cat. #	Description	Size
16672	Tissue-Tack	100ml 500ml

To Order:

In The U.S. Call: 1-800-523-2575 • 215-343-6484
In The U.S. FAX: 1-800-343-3291 • 215-343-0214
In Germany Call: (49) 6221-765767
In Germany FAX: (49) 6221-764620

Should any of our materials fail to perform to our specifications, we will be pleased to provide replacements or return the purchase price. We solicit your inquiries concerning all needs for life sciences work. The information given in this bulletin is to the best of our knowledge accurate, but no warranty is expressed or implied. It is the user's responsibility to determine the suitability for his own use of the products described herein, and since conditions of use are beyond our control, we disclaim all liability with respect to the use of any material supplied by us. Nothing contained herein shall be construed as a recommendation to use any product or to practice any process in violation of any law or any government regulation.