

BioMag[®] SelectaPure[™] Anti-CD4 Antibody

Description

The study of CD4 mediated events is an area of active research in regards to their role in autoimmune disorders, most notably those relating to Human Immunodeficiency Virus (HIV) as CD4 is a receptor for the Human Immunodeficiency Virus type I (HIV-1) envelope protein gp120. The CD4 antigen is expressed on a subset of peripheral blood lymphocytes, the T "helper" lymphocytes. It is also present on most thymocytes and is frequently co-expressed with CD8. CD4 is also expressed on all monocytes, but at a lower density than on CD4+ T lymphocytes. CD4+ T lymphocytes are active in inducing and helping the synthesis of immunoglobulins by B cells. Autoimmune disorders have been shown to be associated with T cell subset abnormalities due to a loss of inducer cells or due to the presence of abnormal T cells.

The CD4 antigen is a single-chain transmembrane glycoprotein with a 59 kDa molecular weight. It binds to a non polymorphic region of MHC Class II molecules. CD4 is a co-receptor in MHC Class II restricted antigen-induced activation.

BioMag Anti-CD4 particles are designed for positive selection of T helper lymphocytes and monocytes. The Anti-CD4 MT310 clone used with BioMag recognizes human T helper lymphocytes and monocytes. BioMag Anti-CD4 particles may be used in the detection of T cell acute lymphoblastic leukemia and lymphoblastic lymphoma and in the identification of T cell chronic lymphocytic leukemia derived from a helper population. BioMag Anti-CD4 is a suspension of magnetic particles approximately 1 µm in size that are attached to a mouse anti-human T cell IgG1 monoclonal antibody.

Figure A





Figure B



General Recommendation*:

Conc. #	1.40 x 10 ⁸ bead/ml
Vol. Used	0.018 ml
# Particles	2.52 x 10 ⁶ per test
# Target Cells	6.77 x 10 ⁵ per test
Bead:Cell Ratio	3.7
% Depletion	99.44%

*These values should be used as a starting point in optimizing experimental protocols. Due to differences in the distribution of cell types in samples and other variables, the researcher is strongly encouraged to determine the optimal particle to cell ratios for their experiments.

Cell sorting results using BioMag Selectapure CD4 anti-human leukocyte particles for positive selection. Typically, whole blood or purified leukocytes and particles are incubated for 30 minutes at room temperature and then magnetically separated. The supernatant is collected, incubated with the appropriate two-color antibody cocktail, and then analyzed by flow cytometry. Figure A depicts the cell population prior to positive selection. Figure B depicts the cell population after positive selection. The particle to cell ratios reported above are based on experiments where cells were exposed to the particles once.

Particle Concentration

The concentration of BioMag is 1.5 mg/ml. There are approximately 1 x 10⁸ BioMag particles per mg.

Cell Separation Recommendations

Depending upon antigen availability and the size of the target cell population, cell sorting applications may require up to 50-60 magnetic particles per cell based on the target cell population. Magnetic particles and cells should be incubated at room temperature for 30 minutes to one hour in media containing 5-10% protein (to reduce non-specific binding) for successful separation. Gentle end over end or rocking during incubation is required for optimal results. (Note: Increasing the incubation time beyond one hour may be necessary to achieve the desired depletion.) Each researcher must optimize particle to cell ratio and incubation time for the application.

Some applications require the detachment of BioMag antibody particles from cells after separation. One approach would involve culturing cells after positive selection. Cultures can be maintained for about 48 hours during which magnetic particles fall away from cells due to cell surface changeover. The magnetic particles are then easily removed via a magnetic separation. Another approach is the use of a protease such as chymopapain to break the antigen-antibody bond and remove the particles magnetically. Depending upon the application, it may not be necessary to remove the cells from the BioMag particles. BioMag particles are only 1µm in size and have been successfully used in FACS equipment. They will not jam the machine and are distinguishable from cells. Alternatively, negative selection approaches can be very effective in producing specific cell populations.

Storage and Stability

The suspension is supplied in PBS/EDTA/1.0% BSA/0.1% sodium azide buffer at pH 7.5. Washing BioMag Anti-CD4 particles in sterile media to remove preservative prior to use is recommended. Using a magnetic separation unit for washing instead of centrifugation is strongly recommended. Do not freeze, dry or centrifuge BioMag. Freezing, drying and centrifuging BioMag can result in aggregation and loss of binding activity. BioMag Anti-CD4 particles are stable when stored at 4°C.

Safety

BioMag Anti-CD4 particle suspension contains sodium azide. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Upon disposal of material, flush with a large volume of water to prevent azide accumulation. Please consult the Material Safety Data Sheet for more information.

Ordering Information:

Cat. #	Description	Size
85004	BioMag Anti-CD4 Antibody	1 ml
		5 ml

To Order:

In The U.S. Call: 1-800-523-2575 • 215-343-6484

In The U.S. FAX: 1-800-343-3291 • 215-343-0214

In Germany Call: (49) 6221-765767

In Germany FAX: (49) 6221-764620

Order online anytime at www.polysciences.com

This product is for research use only and is not intended for use in humans or for in vitro diagnostic use.

Should any of our materials fail to perform to our specifications, we will be pleased to provide replacements or return the purchase price. We solicit your inquiries concerning all needs for life sciences work. The information given in this bulletin is to the best of our knowledge accurate, but no warranty is expressed or implied. It is the user's responsibility to determine the suitability for his own use of the products described herein, and since conditions of use are beyond our control, we disclaim all liability with respect to the use of any material supplied by us. Nothing contained herein shall be construed as a recommendation to use any product or to practice any process in violation of any law or any government regulation.