

TECHNICAL DATA SHEET 602

Page 1 of 2

Periodic Acid Schiff's Stain Kit (PAS)

Catalog Number 24200

Please refer to the MSDS for chemical and safety information.

Introduction

PAS techniques are used to demonstrate polysaccharides, neutral mucosubstances and basement membranes primarily in tissue. The PAS reagent is also called Fielgen Stain for the demonstration of DNA with a different protocol. Kidney is the most sensitive control. The demonstration of glycogen is best represented by a section of liver with a digestion step used as a negative control in the staining.

The principle of the PAS reaction is the conversion or loss of the quinoid structure and the masking of chromophores. This forms a colorless compound called leuco-fuchsin which can be changed to the classic red reaction by washing in running water to remove the sulphurous groups and restore the quinoid groups. Excess Schiff reagent is removed by potassium metabisulfite rinses. Thus preventing a false positive by oxidation of the reagent in the tissue.

Fixation

10% Neutral Buffered Formalin or alcohol can be used. Sections should be cut at 4 μ to 5 μ .

Deparaffinize

Slides should be deparaffinized through xylene or xylene substitute to remove the paraffin from the section and descending grades of alcohol to distilled water just prior to staining.

Staining Procedure: Routine Procedure for Room Temperature

1. Place slides in 0.5% Periodic Acid for 5 minutes
2. Wash slides in 3 changes of distilled water.
3. Place slides in Schiff's Reagent for 15 minutes. The Schiff's Reagent should be brought to room temperature prior to staining. Clean staining dishes or Coplin jars are required. The glassware should be acid cleaned to avoid contamination.
4. Wash for 1 minute each in 2 changes of 0.55% Potassium Metabisulfite to remove excess reagent.
5. Wash in running tap water for 10 minutes to allow the color to develop.
6. Counterstain with Acidified Harris Hematoxylin for 30 seconds.
7. Wash the slides in running tap water to remove excess hematoxylin.
8. Dehydrate through 95% alcohol, and 100% alcohol to xylene.
9. Coverslip with Poly-Mount® or other media.

Staining Procedure: Microwave Procedure

Steps 1 to 4 must be done under a hood and the microwave must be vented or under the hood. Larger amounts of stain can be made and heated by extending the times as indicated in the procedure. Heating in a Coplin jar will require 15 to 30 seconds to reach 60°C and 1 to 1 1/2 minutes for 250 ml to reach 60°C in a larger staining dish on HIGH in a microwave or as directed by the manufacturer using a probe. Please calibrate your oven by checking this temperature with room temperature distilled water for accurate control.

1. The glassware should be acid cleaned to avoid contamination.
2. Preheat the 0.5% Periodic Acid Solution in the microwave for 15 to 30 seconds for 50ml in a Coplin jar or 1 to 1 1/2 minutes for 250mL of solution in a large staining dish. Place slides in warm solution for 2 minutes.
3. Wash slides in 3 changes of distilled water.
4. Preheat the Schiff's reagent for 15 to 30 seconds with 50 ml in a Coplin jar or 1 to 1 1/2 minutes with 250 ml in a large staining dish prior to staining. Place the slides in the Schiff's reagent for 1 to 2 minutes.
5. Wash for 1 minute each in 2 changes of 0.55% Potassium Metabisulfite to remove excess reagent.
6. Wash in running tap water for 10 minutes to allow the color to develop.
7. Counterstain with Acidified Harris Hematoxylin for 30 seconds.
8. Wash the slides in running tap water to remove excess hematoxylin.
9. Dehydrate through 95% alcohol, and 100% alcohol to xylene.
10. Coverslip with Poly-Mount® or other media.

PAS with Amylase Digestion

The digestion should be completed and the slides stained with the patient or test slides. The glycogen will not stain after digestion.

1. Use 100mg of Amylase, provided, to 50ml of pH 5.0 Phosphate Citrate Buffer. Stir until dissolved. Use immediately.
2. Place slides in the Amylase solution for 20 minutes. This solution can be warmed in the microwave however; it will change the pH and must be watched carefully to avoid over digestion of the tissue.
3. Wash in running tap water for 1 minute.
4. Stain as above for negative control.

Results: Glycogen - Red Nuclei - Blue

Ordering Information:

Cat. #	Description	Size
24200-1	Periodic Acid Schiff's Stain Kit (PAS) Kit Contains:	
	0.5% Periodic Acid Aq	250ml
	Schiff' Reagent	250ml
	0.55% Potassium Metabisulfite	250ml
	Harris Hematoxylin Acidified	250ml
	Phosphate Citrate Buffer pH 5.0	250ml
	Diastase Powder	.5gm

Also Available:

08381-120	Poly-Mount®	120ml
08381-940	Poly-Mount®	940ml

To Order:

In The U.S. Call: 1-800-523-2575 • 215-343-6484
 In The U.S. FAX: 1-800-343-3291 • 215-343-0214
 In Germany Call: (49) 6221-765767
 In Germany FAX: (49) 6221-764620

Order online anytime at www.polysciences.com

Should any of our materials fail to perform to our specifications, we will be pleased to provide replacements or return the purchase price. We solicit your inquiries concerning all needs for life sciences work. The information given in this bulletin is to the best of our knowledge accurate, but no warranty is expressed or implied. It is the user's responsibility to determine the suitability for his own use of the products described herein, and since conditions of use are beyond our control, we disclaim all liability with respect to the use of any material supplied by us. Nothing contained herein shall be construed as a recommendation to use any product or to practice any process in violation of any law or any government regulation.