

U.S. Corporate Headquarters 400 Valley Rd. Warrington, PA 18976 1(800) 523-2575 / (215) 343-6484 1(800)343-3291 fax info@polysciences.com Polysciences Europe GmbH Badener Str. 13 69493 Hirschberg an der Bergstrasse, Germany +(49) 06201-845200 +(49) 06201-8452020 fax info@nolysciences de Polysciences Asia-Pacific, Inc. 2F-1, 207 DunHua N. Rd. Taipei, Taiwan 10595 (886) 2 8712 0600 (886) 2 8712 2677 fax info@polysciences.tw

# **TECHNICAL DATA SHEET 559**

Page 1 of 1

# BioMag® Goat anti-Rabbit IgG

Catalog Number: 84300

#### DESCRIPTION

BioMag® Goat anti-Rabbit IgG is a suspension of BioMag® particles approximately 1.5µm in size, which are covalently coated with goat anti-Rabbit IgG antibody. The suspension is supplied in phosphate buffered saline (pH 7.4) with EDTA and sodium azide added. After shaking vigorously or vortexing, BioMag® Goat anti-Rabbit IgG is ready to use.

### **CHARACTERISTICS**

 $\begin{array}{ll} \mbox{Mean Diameter:} & \sim 1.5 \mu\mbox{m} \\ \mbox{Particle Concentration:} & 1 \mbox{ mg/ml} \end{array}$ 

Binding Capacity: 1ml (1mg) of BioMag® Goat anti-Rabbit IgG will

bind > 0.1mg of rabbit lgG.

## **PROCEDURE**

Researchers are advised to optimize the use of BioMag® in any application as procedures designed by other manufacturers may not be ideal.

BioMag® Goat anti-Rabbit IgG is suitable for use as a secondary antibody in radioimmunoassays and enzyme immunoassays that utilize a rabbit primary antibody. In general, BioMag® Goat anti-Rabbit IgG may be used to quickly and conveniently separate a rabbit antibody or rabbit antibody complex from solution. (For additional information on the use of this product in cell separation, please see Technical Data Sheet 528, BioMag® and Cell Sorting.) Separation may be accomplished magnetically or through centrifugation.

#### STORAGE AND SAFETY

Storage Store at 4°C. Freezing, drying, or centrifuging BioMag® may

result in irreversible aggregation and loss of binding activity.

**Safety** This particle suspension contains sodium azide. Sodium azide may react with lead and copper plumbing to form explosive

may react with lead and copper plumbing to form explosive metal azides. Upon disposal of material, flush with a large volume of water to prevent azide accumulation. Please consult

the Safety Data Sheet for more information.

This product is for research use only and is not intended for use in humans or for *in vitro* diagnostic use.

## **ORDERING INFORMATION**

Cat. #	Description	Size
84300-50	BioMag® Goat anti-Rabbit IgG	50ml
84300-500	BioMag® Goat anti-Rabbit IgG	500ml

## **TO ORDER**

In The U.S. Call: 1(800) 523-2575 • (215) 343-6484 In The U.S. Fax: 1(800) 343-3291 • (215) 343-0214

In Germany Call: +(49) 06201-845200 In Germany Fax: +(49) 06201-8452020 In Asia Call: (886) 2 8712 0600 In Asia Fax: (886) 2 8712 2677

Order online anytime at www.polysciences.com