

PSVue™ 380 Reagent Kit

PSVue™ 380 REAGENT KIT CONTENTS

Vial Containing preweighted amount of apo-PSS380 solid dye (at least 0.5 mg)
Vial of 8.4 mM zinc Nitrate solution in water (0.5 mL)

Note: Absolute ethanol is needed to formulate the product but not provided

DESCRIPTION

The **PSVue™ 380** (formerly PSS-380) reagent kit contains components to provide a 1 mM solution of **PSVue™ 380** in aqueous solution. The structure and spectral properties of **PSVue™ 380** are shown in Figures 1 and 2 respectively. **PSVue™ 380** was originally developed as a sensor for phosphate groups [1] but has also been found to bind strongly to phosphatidylserine (PS) residues exposed on membrane surfaces, through its zinc(II)-dipicolylamine (Zn-DPA) functionality, thus making it a useful reagent for detection of apoptosis in different cell lines [2, 3, 4, 5], as well as for monitoring PS levels in other cellular processes [6]. Staining selectivity is very similar to that obtained with annexin V-FITC since the Zn-DPA PS affinity group has been shown to bind to the same membrane sites as annexin V [2]. Microscopy data indicate that **PSVue™ 380** staining is more intense than annexin V, which is attributed to annexin V's much larger size per fluorophore [2]. Negatively charged bacterial cell walls of both Gram-positive and Gram-negative bacteria (e.g. *S. aureus*, *E. coli*) have also been shown to be labeled selectively with **PSVue™ 380** [7] even in the presence of mammalian cells. Given the ubiquitous presence of anionic membranes in bacteria, **PSVue™ 380** is expected to label most strains. An attractive feature of **PSVue™ 380** is that its fluorescence emission increases by almost 10-fold upon binding to a bilayer membrane [7] thus eliminating the need to wash the cells after addition of the fluorophore. In addition to its utility in cell biology research, **PSVue™ 380** may be useful in the automation of biotechnology processes and high-throughput screening

FIGURE 1

PSVue™ 380 and Precursor apo-PSS380

PSVue™ 380 Chemical Data

Molecular Formula: $C_{40}H_{36}N_{10}O_{12}Zn_2$
Molecular Weight: 979.6 g/mol
Extinction Coefficient: $1.69 \times 10^4 M^{-1} cm^{-1}$ (in water)
Quantum Yield: 0.50 (in water)





FIGURE 2

PSVue™ 380 Absorption and Fluorescence Emission Spectra (5 μM in water; abs. max=380 nm; fl.em max=440 nm).



STORAGE/STABILITY

- For long term storage, the kit maybe refrigerated at 4-8°C. Bring to room temperature before use.
- Once formulated the **PSVue™ 380** dye stock must be protected from bright direct light and examined for crystals prior to use. If crystals are noted in the dye stock, it can be warmed slightly to 40°C in a water bath and sonicated or vortexed to redissolve the crystals.
- The **PSVue™ 380** 2mM stock solution should be stored at 4°C and is best used within 5 days.

PROCEDURE**Formulation Procedure to Prepare 2mM stock solution:**

1. Using pre-weighed apo-PSS-380 solid supplied, prepare a 4mM solution of apo-PSS-380 (i.e. 2.404 mg/mL) in absolute ethanol in the 2mL vial. [Note: Make sure the solid is fully dissolved. Sonication for about 20 minutes is necessary].
2. Add an equal volume of the 8.4mM zinc nitrate solution in water provided to the apo-PSS-380 solution from step 1.
3. Sonicate and/or heat the solution in a water bath at 40°C for 30 minutes to provide a homogeneous solution.
4. Label the vial as 2mM **PSVue™ 380** stock solution in 1:1 ethanol/water.

In Vitro Cell Staining Conditions:

1. Typical concentrations of **PSVue™ 380** used for in vitro cell labeling studies are in the range of 1-50 µM.
2. The recommended buffer for cell staining is a TES [N-tris-(hydroxymethyl)-methyl-2-aminoethane sulfonic acid] buffer system comprising (5 mM TES, 145 mM NaCl, pH=7.4), as used in references [8], [9] and [10]. TES buffer should also be used for any wash steps after labeling, if necessary.

Notes:

- i. PBS buffer can cause problems with in vitro cell staining using PSVue dyes due to the presence of anionic phosphate therefore it should NOT be used.
- ii. PSVue 380 binding to apoptotic cells is almost instantaneous [3]

In Vivo Imaging Conditions:

Fluorescence microscopy of **PSVue™ 380** label cells can be performed on a fluorescent microscope with excitation at 350nm and a DAPI/Hoechst/AMCA filter set [2]. Flow cytometry of labeled cells has been performed with an Enterprise II laser with excitation at 350nm and emission collected at 440nm [2].

REFERENCES

1. Ojida A, Mito-oka Y, Inoue M-A and Hamachi J. First artificial receptors and chemosensors towards phosphorylated peptide in aqueous solution. *J. Am. Chem. Soc.*, 2002, 124, 6256-58
2. Koulov AV, Hanshaw RG, Stucker KA, Lakshmi C and Smith BD. Biophysical studies of a synthetic mimic of the apoptosis-detecting protein annexin V. *Israel. J of Chem.*, 2005, 45, 373-379.
3. Koulov AV, Stucker KA, Lakshmi C, Robinson JP and Smith BD. Detection of apoptotic cells using a synthetic fluorescent sensor for membrane surfaces that contain phosphatidylserine. *Cell Death and Differentiation*, 2003, 10, 1357-1359
4. Manaka J, Kuraishi T, Shiratsuchi A, Nakai Y, Higashida H, Henson P and Nakanishi Y. Draper-mediated and phosphatidylserine-independent phagocytosis of apoptotic cells by drosophila hemocytes/macrophages. *J. Biol. Chem.*, 2004, 279, 46, 48466-48476
5. Hanshaw RG and Smith BD. New reagents for phosphatidylserine recognition and detection of apoptosis. *Bioorg. & Med. Chem.* 2005, 13, 5035-5042
6. Fratti RA, Jun Y, Merz AJ, Margolis N, Wickner W. Interdependent assembly of specific regulatory lipids and membrane fusion proteins into the vertex ring domain of docked vacuoles. *J. Cell Biol.*, 2004, 167, 6, 1087-1098.
7. Leevy, W. M.; Johnson, J. R.; Lakshmi, C.; Morris, J.; Marquez, M.; Smith, B. D. Selective recognition of bacterial membranes by zinc(II)-coordination complexes. *Chem. Commun.* 2006, 1595-1597.
8. DiVittorio KM, Johnson JR, Johansson E, Reynolds AJ, Jolliffe KA and Smith BD. Synthetic peptides with selective affinity for apoptotic cells. *Org. Biomol. Chem.* 2006, 4, 1966-1976.
9. Hanshaw RG, Lakshmi C, Lambert TN, Johnson JR and Smith BD. Fluorescent detection of apoptotic cells by using zinc coordination complexes with a selective affinity for membrane surfaces enriched with phosphatidylserine. *ChemBioChem*, 2005, 6, 2214-2220.
10. Leevy WM, Gammon ST, Jiang H, Johnson JR, Maxwell DJ, Jackson EN, Marquez M, Piwinica-Worms D and Smith BD. Optical imaging of bacterial infection in living mice using a fluorescent near-infrared molecular probe. *JACS*, 2006, 128 (51), 16476-16477.

PSVue™ is a trademark of Molecular Targeting Technologies, Inc. **PSVue™** products are sold under an exclusive license from the University of Notre Dame. US Patent # 7,179,616 and others pending.

This product is offered for research purposes only and is not intended for human therapeutic or diagnostic use.

ORDERING INFORMATION

Cat. #	Description	Sizes
25102	PSVue™ 380 Reagent Kit	1 kit

Order online anytime at www.polysciences.com

TO ORDER

In The U.S. Call:	1(800) 523-2575 • (215) 343-6484
In The U.S. Fax:	1(800) 343-3291 • (215) 343-0214
In Germany Call:	+(49) 6221-765767
In Germany Fax:	+(49) 6221-764620