

U.S. Corporate Headquarters  
400 Valley Rd.  
Warrington, PA 18976  
1(800) 523-2575 / (215) 343-6484  
1(800)343-3291 fax  
info@polysciences.com

Polysciences Europe GmbH  
Badener Str. 13  
69493 Hirschberg an der Bergstrasse,  
Germany  
+(49) 06201-845200  
+(49) 06201-8452020 fax  
info@polysciences.de

Polysciences Asia-Pacific, Inc.  
2F-1, 207 DunHua N. Rd.  
Taipei, Taiwan 10595  
(886) 2 8712 0600  
(886) 2 8712 2677 fax  
info@polysciences.tw

## TECHNICAL DATA SHEET 734

Page 1 of 2

# ViaCheck™ Concentration Controls

Catalog Numbers: 26409, 24627, 24628, 24629 or

SingleShots™ BLIVC50NSS, BLI60NS, BLI70NS, BLI80NS

ATTENTION: PLEASE NOTE THAT OUR RECOMMENDED SETTINGS HAVE CHANGED.

### INTRODUCTION

Analytical instruments such as viability analyzers require a program of routine maintenance and QC to ensure that each instrument yields accurate and consistent results, and that comparable performance is achieved between instruments. ViaCheck™ Viability and Concentration Standards may be used to validate image-based viability instruments before they're commissioned, and to ensure optimum performance on an ongoing basis. The microsphere standards are pre-stained, and need only be loaded into the analyzer for confirmation of live / dead ratios and counts. Non-biological surrogates remove the need for sample preparation, and offer exceptional stability and reproducibility.

### DESCRIPTION

ViaCheck Viability and Concentration Controls are comprised of individual or mixed populations of (viable) and (non-viable) microsphere surrogates. The suspension simulates a sample of live / dead (or dying) cells stained with Trypan Blue. ViaCheck are intended to serve as reference materials for use with image-based instruments that rely on the Trypan Blue dye exclusion method. ViaCheck have been successfully used with instruments such as the Vi-CELL®, CEDEX® HiRes, Countess™, etc.

### CHARACTERISTICS

26409 or BLIVC50NSS 0.5 x 10<sup>6</sup> particles / mL  
24627 or BLI60NS 1 x 10<sup>6</sup> particles / mL  
24628 or BLI70NS 4 x 10<sup>6</sup> particles / mL  
24629 or BLI80NS 8 x 10<sup>6</sup> particles / mL

\* Nominal values. Specific value reported on Certificate of Analysis.

### MATERIAL

#### Material Supplied

- Microspheres suspended in a solution of buffered salts and surfactant containing 0.08% sodium azide in a 20mL bottle or SingleShot™ vial.

#### Material Required

- Cell viability analyzer (ex. Coulter Vi-CELL™ XR Cell Viability Analyzer)
- Precision pipets with disposable tips to deliver 500µl
- Isotonic buffered saline diluent (optional)

### PROCEDURE

Researchers are advised to optimize the use of particles in any application. For the best accuracy, be sure to work carefully and quickly when sampling and pipetting ViaCheck™ particles. Allowing the particles to stand for even a short period of time could lead to inaccurate data and results. We suggest that only a single, but no more than 3, samples be loaded into the carousel at a time to safeguard against settling.

1. Gently mix (manual inversion or tube rotator) the vial of particles to ensure a well mixed suspension.
2. Place a 500µL sample of the particles into an analyzer sample cup.
3. Place the sample cup in the analyzer sampling station.
4. Using the Vi-CELL™ XR analyzer menu, set up and save a "CELL TYPE" for Viability controls at the settings below. *Note:* The settings below allow the user to analyze the ViaCheck™ Viability Control Particles and may have to be adjusted for each instrument.

Cell Type	Concentration Settings
Minimum Cell Diameter (µm)	5
Maximum Cell Diameter (µm)	50
Minimum Circularity	0
Dilution Factor	1.0
Cell Brightness (%)	85
Cell Sharpness (%)	100
Viable Cell Spot Brightness (%)	75
Viable Cell Spot Area (%)	5
Decluster Degree	Low
Aspirate Cycles	1
Trypan Blue Mixes	3


5. Analyze the sample according to the analyzer's instructions.

### IMPORTANT NOTE ON EXPECTED RESULTS

Certificates of Analysis (COAs) for ViaCheck™ products provide formal lot-specific values for concentration and viability that may be used to establish instrument QC pass / fail criteria. Facilities may establish specific pass / fail criteria after taking historical instrument performance and study objectives into consideration. Users will often base these criteria around the lot-specific result that is issued, e.g. results within some percentage of the reported values on the COA based on established control limits. See "Establishing control limits for the instrument" in *TSD 0711 ViaCheck™ for Cell Viability Analyzers: Best Practices* on bangslabs.com.

### REFERENCES

Lew C, Gomez JA, Rhyner MN. "Instrument-to-instrument Variability in the Vi-CELL Automated Viability Analyzer." [www.particle.com](http://www.particle.com), 2012, Beckman Coulter



**This product is for research use only and is not intended for use in humans or for *in vitro* diagnostic use.**

**ORDERING INFORMATION**

Cat. #	Description	Size
24609	ViaCheck™ Concentration Control (0.5 x 10 <sup>6</sup> )	20ml
BLIVC50NSS	ViaCheck™ Concentration Control (0.5 x 10 <sup>6</sup> ) SingleShots™	25 or 75 vials
24627	ViaCheck™ Concentration Control (1 x 10 <sup>6</sup> )	20ml
BLI60NS	ViaCheck™ Concentration Control (1 x 10 <sup>6</sup> ) SingleShots™	25 or 75 vials
24628	ViaCheck™ Concentration Control (4 x 10 <sup>6</sup> )	20ml
BLI70NS	ViaCheck™ Concentration Control (4 x 10 <sup>6</sup> ) SingleShots™	25 or 75 vials
24629	ViaCheck™ Concentration Control (8 x 10 <sup>6</sup> )	20ml
BLI80NS	ViaCheck™ Concentration Control (8 x 10 <sup>6</sup> ) SingleShots™	25 or 75 vials

**RELATED PRODUCTS**

Cat. #	Description	Size
24622	ViaCheck™ 0% Viability Control	20ml
BLI10BS	ViaCheck™ 0% Viability Control SingleShots™	25 or 75 vials
25997	ViaCheck™ 25% Viability Control	20ml
BLI25BS	ViaCheck™ 25% Viability Control SingleShots™	25 or 75 vials
24623	ViaCheck™ 50% Viability Control	20ml
BLI20BS	ViaCheck™ 50% Viability Control SingleShots™	25 or 75 vials
24624	ViaCheck™ 75% Viability Control	20ml
BLI30BS	ViaCheck™ 75% Viability Control SingleShots™	25 or 75 vials
24625	ViaCheck™ 90% Viability Control	20ml
BLI40BS	ViaCheck™ 90% Viability Control SingleShots™	25 or 75 vials
24626	ViaCheck™ 100% Viability Control	20ml
BLI50BS	ViaCheck™ 100% Viability Control SingleShots™	25 or 75 vials

RESULTS	
Cell Count	1156
Viable Cell Count	0
Viability (%)	0
Total Cells / mL (x 1.0E6)	1.03
Viable Cells / mL (x 1.0E6)	0
Average Diameter (µm)	10.25
Average Circularity	0.95
Images	50
Average Cells / Image	23.1
Average Background Intensity	205

Representative Vi-CELL® XR data : ViaCheck™ Concentration Control Particles

**RELATED LITERATURE**

TSD 0711 ViaCheck™ for Cell Viability Analyzers: Best Practices  
 TSD 1008 Handling & Pipetting Concentration Standards  
 TSD 0708 Optimization of Vi-CELL® XR Settings for calibration using ViaCheck™ Controls.

**TRADEMARKS AND REGISTERED TRADEMARKS**

ViaCheck™ & SingleShots™ are trademarks of Polysciences, Inc.  
 Vi-CELL® is a registered trademark of Beckman Coulter, Inc.  
 CEDEX® is a registered trademark of Roche Inc.  
 Countess™ is a trademark of Thermo Fisher.

**STORAGE AND SAFETY**

**Storage** Store at 2-8°C. Refrigerated storage is intended to deter the growth of opportunistic microorganisms within the suspensions; it is important to note that biocontamination would alter particle counts. Freezing of particles may result in irreversible aggregation.

**Safety** This particle suspension contains sodium azide. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Upon disposal of material, flush with a large volume of water to prevent azide accumulation. Please consult the Safety Data Sheet for more information.

**TO ORDER**

In The U.S. Call: 1(800) 523-2575 • (215) 343-6484  
 In The U.S. Fax: 1(800) 343-3291 • (215) 343-0214

In Germany Call: +(49) 06201-845200  
 In Germany Fax: +(49) 06201-8452020

In Asia Call: (886) 2 8712 0600  
 In Asia Fax: (886) 2 8712 2677

Order online anytime at [www.polysciences.com](http://www.polysciences.com)

Should any of our materials fail to perform to our specifications, we will be pleased to provide replacements or return the purchase price. We solicit your inquiries concerning all needs for life sciences work. The information given in this bulletin is to the best of our knowledge accurate, but no warranty is expressed or implied. It is the user's responsibility to determine the suitability for their own use of the products described herein, and since conditions of use are beyond our control, we disclaim all liability with respect to the use of any material supplied by us. Nothing contained herein shall be construed as a recommendation to use any product or to practice any process in violation of any law or any government regulation.