

U.S. Corporate Headquarters 400 Valley Rd. Warrington, PA 18976 1(800) 523-2575 / (215) 343-6484 1(800)343-3291 fax info@polysciences.com Polysciences Europe GmbH Handelsstrasse 3 D-69214 Eppelheim, Germany +(49) 6221-765767 +(49) 6221-764620 fax info@polysciences.de Polysciences Asia-Pacific, Inc. 2F-1, 207 Dunhua N. Rd. Taipei, Taiwan 10595 (886) 2 8712 0600 (886) 2 8712 2677 fax info@polysciences.tw

## **TECHNICAL DATA SHEET 259**

# **Macroinvertabrate Mounting Medium**

Page 1 of 1

CMCP is a non-resinous, water miscible mounting media used to make permanent transparent mounts. With adequate clearing of the specimens, this medium can be used to mount almost any material, provided it is not stained, calciferous material, or thin tissue sections. Although used like Canada balsam, this embedding media has the advantage of allowing live or preserved specimens to be mounted directly from water or alcohol.

Free living nematodes and live parasitic worms are among the organisms suitable for mounting with CMCP. The CMCP mounting media provide them with rapid relaxation, death, and clearing action.

Whereas specimens with little air can be mounted alive, when they possess extensive air passages and/or large air sacs, immersion in a low viscosity medium is necessary to displace trapped air. To avoid trapping air when CMCP-9 macroinvertebrate mounting medium is used, allow the medium to evaporate for a 24 hour period before attempting to coverslip. Specimens made with the CMCP-10 should dry for one week before adding a drop of CMCP-10 and coverslipping.

If the specimen is thick, it may be desirable to circle the coverslip with the mounting media. Several rings of mounting media may be needed at day intervals to achieve the proper height for coverslipping. If staining is desirable after mounting, this can be completed within 24 hours.

#### **WARNING:**

The full chemical, physical, and toxicological properties for this product have not been established. Some components cause irritation and may also be skin sensitizers.

The presence of phenol makes this mixture toxic and irritating. Phenol is a skin, eye, and respiratory system irritant; its solution can be absorbed rapidly through the skin and death can result from collapse within 30 minutes to several hours. Spills on the skin can cause headaches, dizziness, muscular weakness, dimness of vision, ringing in the ears, irregular and rapid breathing, weak pulse, and dyspnea (difficult/labored breathing) may also develop, and may be followed by loss of consicousness, collapse, and death. Lactic acid is also a skin irritant.

### **ORDERING INFORMATION:**

Cat. # Description Size

16299 CMCP-9 Mounting Media Low viscosity, colorless mountant 100mL, 500mL

Refractive index about 1.40

16300 CMCP-10 Mounting Media High Viscosity, colorless mountant 250mL, 500mL

Refractive index about 1.41

Some researchers like to tint their mounting media with Aniline blue.

O2570 Aniline blue (water soluble), C.I. 42755 certified 25g

#### TO ORDER

In The U.S. Call: 1(800) 523-2575 • (215) 343-6484 In The U.S. Fax: 1(800) 343-3291 • (215) 343-0214

In Germany Call: +(49) 6221-765767 In Germany Fax: +(49) 6221-764620

In Asia Call: (886) 2 8712 0600 In Asia Fax: (886) 2 8712 2677

Order online anytime at www.polysciences.com