



U.S. Corporate Headquarters
400 Valley Rd.
Warrington, PA 18976
1(800) 523-2575 / (215) 343-6484
1(800)343-3291 fax
info@polysciences.com

Polysciences Europe GmbH
Badener Str. 13
69493 Hirschberg an der Bergstrasse,
Germany
+(49) 06201-845200
+(49) 06201-8452020 fax
info@polysciences.de

Polysciences Asia-Pacific, Inc.
2F-1, 207 DunHua N. Rd.
Taipei, Taiwan 10595
(886) 2 8712 0600
(886) 2 8712 2677 fax
info@polysciences.tw

TECHNICAL DATA SHEET 530

Page 1 of 2

BioMag[®] Streptavidin, Nuclease-free

Catalog Number: 8MB4804

DESCRIPTION

BioMag[®] Streptavidin is a nuclease-free suspension of BioMag[®] particles approximately 1.5µm in size, which are covalently coated with streptavidin. The suspension is supplied in a phosphate buffered saline (pH 7.4) containing 0.1% BSA. Sodium azide has been added as an anti-microbial. Shake vigorously or vortex before use. Magnetically separate the BioMag[®] particles, aspirate the supernatant, and resuspend in an appropriate buffer.

CHARACTERISTICS

Mean Diameter: ~1.5µm
Particle Concentration: 1mg/ml
Binding Capacity: 1mg of BioMag[®] Streptavidin will bind:
>1500 pmoles of free biotin
>1000 pmoles of a 20-mer biotinylated oligonucleotide
>200 pmoles of a 100-mer biotinylated oligonucleotide
>70 pmoles of a 300 bp 5-biotinylated double-stranded DNA
>25 pmoles of a 1Kbp 5-biotinylated double-stranded DNA

MATERIAL

Material Supplied

- BioMag[®] Streptavidin: 10ml or 25ml

Material Required

- Binding Buffer: 20mM Tris and 0.5M NaCl at pH 8.0
- Wash Buffer: 7mM Tris and 0.17M NaCl at pH 8.0
- DEPC-treated water
- Nuclease-free microcentrifuge tubes
- Magnetic separator

PROCEDURE

Researchers are advised to optimize the use of BioMag[®] in any application as procedures designed by other manufacturers may not be ideal.

The following procedure is for the isolation of 1-2µg of mRNA from approximately 75-100µg of total RNA. The total isolation time is less than 30 minutes.

1. Dispense 200µl of BioMag[®] Streptavidin into a nuclease-free microcentrifuge tube. Using a magnetic separation unit, pull the magnetic particles to the side of the microcentrifuge tube for 30 seconds. Remove and discard the supernatant. Resuspend the

BioMag[®] Streptavidin in 100µl of Binding Buffer.

2. Incubate 2.5µl (2.5µg) of 5-Biotinylated Oligo (dT) (or an appropriate amount of biotinylated molecule) with the 100µl of BioMag[®] Streptavidin from Step 1 for 15 minutes at room temperature.
3. Magnetically separate for 30 seconds and discard the supernatant. Wash the Oligo (dT) bound particles from Step 2 with 100µl of Binding Buffer 2 times, leaving the magnetic particles as a wet cake.
4. Bring up the total RNA sample in DEPC-treated water to a total volume of 90µl.
5. Incubate the RNA sample at 55°C for 5 minutes to disrupt secondary structures.
6. Add 10µl of 5M NaCl to achieve a final concentration of 0.5M NaCl.
7. Add the total RNA to the washed magnetic particles from Step 3. Mix gently and hybridize at room temperature for 3 minutes.
8. Magnetically separate and wash the particles with 100µl of Wash Buffer.
9. Elute the bound mRNA with 25-50µl of DEPC-treated water at 55°C for 2 minutes.
10. Magnetically separate and transfer the supernatant to a nuclease-free microcentrifuge tube.
11. Repeat elution of mRNA with 25-50µl of DEPC-treated water at 55°C for another 2 minutes in order to completely elute the bound mRNA from the particles. Magnetically separate and transfer the supernatant to the tube containing the first elution of mRNA from Step 10.

REFERENCES

1. **Hornes, E., L. Korsnes.** 1990. Magnetic DNA hybridization properties of oligonucleotide probes attached to superparamagnetic beads and their use in the isolation of poly(A) mRNA from eukaryotic cells. *Genet Anal Tech Appl*; 7(6):145-150.
2. **Morrissey, D.V., M. Lombardo, J.K. Eldredge, K.R. Kearney, E.P. Groody, M.L. Collins.** 1989. Nucleic acid hybridization assays employing dA-tailed capture probes. Multiple capture methods. *Anal Biochem*, 181(2):345-359.

STORAGE AND STABILITY

Store at 4°C. Freezing, drying, or centrifuging BioMag[®] may result in irreversible aggregation and loss of binding activity.

SAFETY

This particle suspension contains sodium azide. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Upon disposal of material, flush with a large volume of water to prevent azide accumulation. Please consult the Safety Data Sheet for more information.

These products are for research use only and are not intended for use in humans or for *in vitro* diagnostic use.

ORDERING INFORMATION

Cat. #	Description	Sizes
8MB4804-10	BioMag® Streptavidin, Nuclease-free	10ml
8MB4804-25	BioMag® Streptavidin, Nuclease-free	25ml
8MB4804-1	BioMag® Streptavidin, Nuclease-free	100ml

TO ORDER

In The U.S. Call: 1(800) 523-2575 • (215) 343-6484

In The U.S. Fax: 1(800) 343-3291 • (215) 343-0214

In Germany Call: +(49) 06201-845200

In Germany Fax: +(49) 06201-8452020

In Asia Call: (886) 2 8712 0600

In Asia Fax: (886) 2 8712 2677

Order online anytime at www.polysciences.com