

Methylene Blue - Basic Fuchsin

Staining Procedure for Plastic Embedded Tissue

Working Solution:

Methylene Blue, Basic Fuchsin

Procedure:

1. Immerse slides in Methylene Blue-Basic Fuchsin at 65°C.
The time required varies depending on the thickness and type of embedding media.
2. Rinse in 2 changes of distilled water.
3. Stain in Basic Fuchsin at room temperature for 20 minutes.
The time required varies depending on the thickness and type of embedding media.
4. Rinse thoroughly in distilled water.
5. Dry stained sections on a 40°C hot plate.
6. Mount sections.

Stain Results:

Nuclei - blue
Cytoplasm, Mitochondria, Cilia - reds or pinks
Cartilage / Cell granules - blue to purple

Product Information:

Cat. #	Size	Description
04220	50g	Methylene blue chloride, C.I. 52015, certified
06342	25g	Fuchsin, basic, C.I. 42500, certified

To Order:

In The U.S. Call: 1-800-523-2575 • 215-343-6484
In The U.S. FAX: 1-800-343-3291 • 215-343-0214

In Germany Call: (49) 6221-765767
In Germany FAX: (49) 6221-764620