



U.S. Corporate Headquarters
400 Valley Rd.
Warrington, PA 18976
1(800) 523-2575 / (215) 343-6484
1(800)343-3291 fax
info@polysciences.com

Polysciences Europe GmbH
Handelsstrasse 3
D-69214 Eppelheim, Germany
+(49) 6221-765767
+(49) 6221-764620 fax
info@polysciences.de

Polysciences Asia-Pacific, Inc.
2F-1, 207 Dunhua N. Rd.
Taipei, Taiwan 10595
(886) 2 8712 0600
(886) 2 8712 2677 fax
info@polysciences.tw

TECHNICAL DATA SHEET 305B

Page 1 of 1

LR White Special Handling Instructions

Sometimes due to warm weather and length of time in shipment, it may be necessary to ship LR White with dry ice to keep the integrity of the material. This will cause freezing but is not harmful to the performance of the resin.

Different monomers will freeze out of the solution at different temperatures, leading to a slush form from about -20°C downwards.

After receiving material, keep the resin at room temperature long enough to allow the material to thaw, then place immediately into a refrigerator at 4°C. When ready to use gently heat material in a water bath set at 35-40°C and mix until the monomers go back into a uniform solution. Prolonged exposure to warm temperatures will cause an increase in viscosity, harming the excellent penetrating properties of this resin.