

Verhoeff Van Gieson Elastin Stain Kit

DESCRIPTION

This method is used for identifying elastic fibers in tissues such as skin, aorta, etc. on formalin-fixed, paraffin-embedded sections, and may be used for frozen sections as well. The elastic fibers will be stained blue-black and background will be stained yellow.

Fixation: 10% formalin.

Section: paraffin sections at 5 µm

KIT CONTAINS

500ml each of the following:

Solution A - 5% ALCOHOLIC HEMATOXYLIN
Solution B - 10% AQUEOUS FERRIC CHLORIDE
Solution C - IODINE SOLUTION

Solution D - 5% SODIUM THIOSULPHATE
Solution E - COUNTERSTAIN I
Solution F - COUNTERSTAIN II

VERHOEFF'S WORKING SOLUTION

The working staining solution should be made up fresh for best results. It will not stain satisfactorily if it is kept more than one working day. Prepare the working solution by adding in order the following reagents:

Solution A ----- 20 ml
Solution B ----- 8 ml
Solution C ----- 8 ml

Mix the above amounts (or needed proportions thereof) well. Solution should be jet black. Use immediately and discard after use.

2% AQUEOUS FERRIC CHLORIDE (PREPARE FRESH)

Solution B ----- 10 ml
Distilled Water ----- 50 ml

PROCEDURE

1. Deparaffinize and hydrate slides to distilled water
2. Stain in Verhoeff's solution for 1 hour, tissue should be black
3. Rinse in tap water with 2-3 changes
4. Differentiate in 2% ferric chloride for 1-2 minutes
5. Stop differentiation with several changes of tap water and check microscopically for black elastic fiber staining and gray background. It is better to slightly under differentiate the tissue, since the subsequent Van Gieson's counterstain can sort of extract the elastic stain
6. Wash slides in tap water
7. Treat with Solution D for 1 minute, discard solution
8. Wash in running tap water for 5 minutes
9. Counterstain in Counterstain I, for 3-5 minutes. Only use counterstain II in place of a counterstain when staining nervous tissue
10. Dehydrate quickly through 95% alcohol, 2 changes of 100% alcohol
11. Clear in 2 changes of xylene for 3 minutes each
12. Coverslip with resinous mounting medium

RESULTS

Elastic fibers ----- blue-black to black
Nuclei ----- blue to black
Collagen ----- red
Other tissue elements ----- yellow

POSITIVE CONTROLS

Aorta, Kidney, Myometrium

SAFETY

Work in a well ventilated area, wear gloves (nitrile), lab coat and goggles. Avoid contact and inhalation. Picric acid: can become explosive if allowed to dry out. Toxic through skin exposure.

STORAGE AND HANDLING

Storage

Store at room temperature

Handling

Glove, chemical goggles & fume hood

ORDERING INFORMATION

Cat. #	Description	Sizes
25089	Verhoeff Van Gieson Elastin Stain Kit	1 kit

Order online anytime at www.polysciences.com

TO ORDER

In the U.S. Call: 1(800) 523-2575 • (215) 343-6484
In the U.S. Fax: 1(800) 343-3291 • (215) 343-0214
Germany Call: +(49) 6221-765767
Germany Fax: +(49) 6221-764620