

Hematoxylin




Our Hematoxylin products offer you the highest quality for routine staining. Strict quality control for all Gold Standard Series Routine Stains ensures lot-to-lot consistency.

Gold Standard Series Gill's Hematoxylin

Gill's Hematoxylin #1 for Cytology

Gill's #1 is ideal for routine cytology staining. This single strength formulation optimally stains gynecological and non-gynecological specimens.



Cheek cells, Gill's #1, 20x

Catalog #	Size
24242-500	500ml
24242-1000	1000ml

Gill's modified EA & Gill's modified OG-6

Orange-6 and EA are the two cytoplasmic counterstain solutions that are used sequentially in the Papanicolaou staining method for clinical cytology. Gill's modified EA and OG-6 are stable in solution and give predictable high quality staining results not previously possible with other formulations.

Gill's modified EA


Catalog #	Size
09783-500	500ml
09783-1000	1000ml
09783-3.75	3.75 liters

Gill's modified OG-6

Catalog #	Size
09782-500	500ml
09782-1000	1000ml
09782-3.75	3.75 liters

Gill's Hematoxylin #2, double strength for Histology & Cytology

Gill's #2 should be used when a stronger or darker nuclear stain is required for cytology samples or for immunohistochemical counterstaining.




Breast, Gill's #2, 10x

Catalog #	Size
24243-500	500ml
24243-1000	1000ml

Harris Hematoxylin, Acidified (mercury-free)

Gold Standard Harris Hematoxylin for routine histology and cytology. Similar results can be achieved to those of Gill's Hematoxylin #1, #2, and #3 formulations by varying staining times.




Colon carcinoma, Harris, 20x

Catalog #	Size
24245-500	500ml
24245-1000	1000ml

Gill's Hematoxylin #3, triple strength for Histology

Gill's #3 is the formulation of choice for routine histological staining when a darker nuclear stain is desired.




Kidney, Gill's #3, 20x

Catalog #	Size
24244-500	500ml
24244-1000	1000ml

Mayers Hematoxylin

For use in both cytology and histology. General purpose nuclear stain, used with hematoxylin and eosin staining. Use as a counterstain for immunohistochemistry procedures.



Colon, Mayers, 20x

Catalog #	Size
24821-500	500ml
24821-1	1 liter