


# 单体产品指南



含酸单体  
丙烯系单体  
粘合促进性单体  
含胺单体  
交联型丙烯酸系单体  
双重活性丙烯酸系单体  
环氧化物/酸酐/亚酰胺  
氟化丙烯酸系单体  
丙烯酸系单体  
高/低折射率单体  
含羟基单体  
单/双官能基乙二醇系寡聚单体  
苯乙烯系单体  
含磺酸盐单体  
紫外光敏感性单体  
乙烯系单体  
用于单体的聚合抑制剂



### Monomers

Polysciences stocks a wide portfolio of monomers. Such variety offers the synthetic chemist the tools to make a rich array of polymer compositions. Our Monomers Selection Guide organizes this set of reactive monomers into various groupings. These organized sets of data will enable the scientist to quickly determine which specific monomer can be used to synthesize custom polymers that meet their needs. In reviewing the data in these selection guides, you will be able to compare and contrast monomer alternatives quickly. More detailed information and chemical structures are included in the alphabetical listing which follows.

### Monomer Selection Guide At-A-Glance

Acid Containing Monomers . . . . .	3
Acid Containing Monomers, Metal Salts . . . . .	3
Acrylic Monomers (Neutral, Monofunctional) . . . . .	4
Adhesion Promoting Monomers . . . . .	5
Amine Containing Monomers . . . . .	6
Crosslinking Acrylic Monomers - Difunctional . . . . .	7
Crosslinking Acrylic Monomers - Multifunctional . . . . .	8
Dual Reactive Acrylic Monomers . . . . .	8
Epoxides / Anhydrides / Imides . . . . .	9
Fluorescent Acrylic Monomers . . . . .	10
Fluorinated Acrylic Monomers . . . . .	11
High / Low Refractive Index Monomers . . . . .	11
Hydroxy Containing Monomers . . . . .	12
Mono and Difunctional Glycol Oligomeric Monomers . . . . .	13
Styrenic Monomers . . . . .	14
Sulfonate Containing Monomers . . . . .	14
UV (light) Active Monomers . . . . .	14
Vinyl and Ethenyl Monomers . . . . .	15
Polymerization Inhibitors for Monomers . . . . .	16

## Acid Containing Monomers

Acidic groups are often used to convey solubility to polymers in aqueous media. These moieties can be converted to a wide range of alternative functional groups. Acid groups can be utilized as catalysts for chemical reactions. Additionally they are employed in polymers as a functional group which enables improved adhesion to a variety of substrates through hydrogen bonding or metal chelation.

	Polymerizable Sites	Polymerization Synthon	Form of Acid	Special Features			
Acrylic acid min. 99.5%	mono	acrylic	carboxylic acid	makes water soluble polymers	00020-250	250 g	128
Acrylic anhydride min. 90%	mono	acrylic	carboxylic acid (protected)	not a crosslinker	00488-50	50 g	128
tert-Butyl methacrylate	mono	acrylic	carboxylic acid (protected)	acid formed by thermal elimination of isobutylene, homopolymer T <sub>g</sub> = 230° C	02058-100	100 g	134
Beta-Carboxyethyl acrylate, >98% Active	mono	acrylic	carboxylic acid	high purity, hydrophilic	24891-100	100 g	131
Methacrylic acid, 99.9%	mono	acrylic	carboxylic acid	high purity, hydrophilic	24897-250	250 g	147
Methacrylic acid, min. 99.5%	mono	acrylic	carboxylic acid	offers latex stability, homopolymer T <sub>g</sub> = 185° C	00212-450	450 g	147
4-Methacryloxyethyl trimellitic anhydride	mono	acrylic	anhydride	adhesion promoter through anhydride	17285-10	10 g	148
Methacryloyl-L-Lysine	mono	acrylic	amino acid	zwitterionic, can derivatize acid or amine	24315-5	5 g	149
o-Nitrobenzyl methacrylate, min. 95%	mono	acrylic	carboxylic acid (protected)	acid formed by photolabile deprotection	24360-10	10 g	151
2-Propene-1-sulfonic acid, sodium salt, 35% soln. in water	mono	vinyl	sulfonic acid salt	water soluble	00064-10	10 g	154
2-Sulfoethyl methacrylate, >90%	mono	acrylic	sulfonic acid ester	water soluble	02597-50	50 g	156
Trichloroacrylic acid	mono	acrylic	carboxylic acid	carboxylic acid	02686-10	10 g	157
4-Vinylbenzoic acid	mono	styrenic	carboxylic acid	aromatic acid	04485-5	5 g	159

## Acid Containing Monomers, Metal Salts

	Polymerizable Sites	Polymerization Synthon	Form of Acid	Special Features			
Barium methacrylate, >95%	dual	acrylic	metal salt	metal salt	01994-50	50 g	130
Lead methacrylate	mono	acrylic	heavy metal salt	x-ray opacity and capture	16382-100	100 g	146
2-ethylhexanoate in methyl methacrylate (0.54/0.46 moles)							
Lithium methacrylate	mono	acrylic	metal salt	water soluble	17117-50	50 g	146-147
Magnesium acrylate	dual	acrylic	metal salt	ionomeric crosslinking, high T <sub>g</sub> polymers	02467-10	10 g	147
Sodium acrylate	mono	acrylic	metal salt	can make high T <sub>g</sub> salt polymers	01207-50	50 g	155
3-Sulfopropyl acrylate, potassium salt	mono	acrylic	sulfonic acid salt	water soluble	17209-100	100 g	156
3-Sulfopropyl methacrylate, potassium salt, 98%	mono	acrylic	sulfonic acid salt	water soluble	17210-100	100 g	156
Zinc (di) methacrylate	dual	acrylic	metal salt	ionomeric crosslinking	03011-100	100 g	160

## Acrylic Monomers (Neutral, Monofunctional)

Monofunctional acrylics shape the type and nature of the main chain polymer backbone. Monomers are chosen to obtain the desired glass transition temperature, flexibility, mechanical strength, polarity and hydrophilic/hydrophobic character of the resulting polymer. Generally, acrylamides exhibit improved resistance to hydrolysis compared to acrylic/methacrylic esters.

	Homopolymer Tg (°C)	Special Features			
Acrylamide, Chemzymes Ultra Pure®	165	hydrophilic, hydrogel synthesis	00019-100	100 g	128
N-Acryloylmorpholine	147	hydrophilic	21192-50	50 g	128
Benzhydryl methacrylate		high refractive index, aromatic	24286-10	10 g	130
Benzyl acrylate, ~99%	6	high refractive index, aromatic	01997-100	100 g	130
Benzyl methacrylate, min. 95%	54	high refractive index, aromatic	02000-100	100 g	130
N-Benzylmethacrylamide		high refractive index, aromatic	17969-25	25 g	130
2-n-Butoxyethyl methacrylate, ~97%		moderate polarity	02034-100	100 g	133
tert-Butyl acrylate	43	hydrocarbon building block monomer	02039-250	250 g	133
sec-Butyl acrylate	12	hydrocarbon building block monomer	02038-25	25 g	133
n-Butyl acrylate, min. 99%	-54	hydrocarbon building block monomer	02037-500	500 g	133
sec-Butyl methacrylate	60	hydrocarbon building block monomer	02057-10	10 g	133
n-Butyl methacrylate	20	hydrocarbon building block monomer	02059-500	500 g	134
iso-Butyl methacrylate	53	hydrocarbon building block monomer	02056-500	500 g	133
tert-Butyl methacrylate	107	hydrophobic / protected acid	02058-100	100 g	134
4-Chlorophenyl acrylate	58	high refractive index, aromatic	01331-10	10 g	134
Cyclohexyl acrylate, min 85%	19	aliphatic, hydrophobic	02109-100	100 g	135
Cyclohexyl methacrylate, ~98%	83	aliphatic, hydrophobic	01837-100	100 g	135
iso-Decyl acrylate	-55	hydrophobic, low Tg	03008-100	100 g	136
n-Decyl methacrylate, 99%	-30	hydrophobic	23344-25	25 g	136
iso-Decyl methacrylate, min. 90%		hydrophobic	22493-100	100 g	136
N,N-Diethylacrylamide, min. 95%		polar building block	00871-25	25 g	137
N,N-Dimethylacrylamide, min. 98%	89	useful for chromatographic resins	02255-100	100 g	138
N,N-Dimethylmethacrylamide		hydrolytic stability, useful for hydrogels	02270-25	25 g	139
N-Diphenylmethylacrylamide		aromatic	19663-5	5 g	139
n-Dodecyl acrylate	-3	hydrophobic, long chain alkyl	02460-50	50 g	140
n-Dodecyl methacrylate	-55	hydrophobic. Low Tg , long chain alkyl	02461-250	250 g	140
N-Dodecylmethacrylamide	15	hydrophobic, hydrolytic stability	04135-10	10 g	140
2-(2-Ethoxyethoxy)ethyl acrylate	-70		02626-100	100 g	140
2-Ethylhexyl acrylate	-50	hydrophobic	00587-250	250 g	141
N-Ethylmethacrylamide		hydrolytic stability, useful for hydrogels	02322-10	10 g	141
1-Hexadecyl methacrylate	22	hydrophobic, aliphatic	02396-25	25 g	143
n-Hexyl acrylate, min. 98%	-45	hydrophobic, aliphatic	02411-100	100 g	144
2-Methoxyethyl acrylate	-50	low Tg monomer	02487-100	100 g	148
2-Methoxyethyl methacrylate, min.85%	16	moderate polarity	02488-100	100 g	148
Methyl Methacrylate, min 99.5%	105	versatile building block monomer	00834-1	1 liter	149
2-Naphthyl acrylate	24	fluorescent, hydrophobic, aromatic	06024-1	1 g	150
N-(n-Octadecyl)acrylamide		hydrophobic	04673-10	10 g	151
n-Octyl methacrylate	-20	hydrophobic	23355-25	25 g	151
N-tert-Octylacrylamide		hydrophobic, hydrolytic stability	03141-25	25 g	152
Neopentyl methacrylate		hydrophobic	06127-10	10 g	152
Pentabromophenyl acrylate		brominated, high RI, flame retardant	06344-10	10 g	152
Pentabromophenyl methacrylate		brominated, high RI, flame retardant	04253-5	5 g	152
Pentafluorophenyl acrylate		perfluorophenyl, low surface energy	06349-5	5 g	152
Pentafluorophenyl methacrylate, 95%		perfluorophenyl, low surface energy	06350-5	5 g	152

			Catalog #	Size	Page
	<b>Homopolymer</b>	<b>Special Features</b>			
	<b>Tg (°C)</b>				
2-Phenoxyethyl methacrylate	54	aromatic, hydrophobic	02640-100	100 g	153
Phenyl acrylate, min. 95%	57	aromatic, UV absorbing	02642-10	10 g	153
Phenyl methacrylate, >95%	110	moderate UV absorbing, aromatic, hydrophobic	02644-10	10 g	153
2-Phenylethyl acrylate, min. 92%	-3	moderate UV absorbing, aromatic, hydrophobic	02834-100	100 g	153
2-Phenylethyl methacrylate, min. 92%	26	moderate UV absorbing, aromatic, hydrophobic	02911-100	100 g	153
n-Propyl acrylate	-37	building block monomer	03132-25	25 g	155
n-Propyl methacrylate	35	building block monomer	03174-100	100 g	155
N-iso-Propylacrylamide	85	hydrophilic	02455-100	100 g	155
Stearyl acrylate	35	hydrophobic, can form crystal domains	02636-100	100 g	155
Stearyl methacrylate	38		02637-100	100 g	156
Tetrahydrofurfuryl acrylate	-12		02907-250	250 g	156
Tribromoneopentyl methacrylate		bromo non-aromatic	03057-10	10 g	157
2,4,6-Tribromophenyl acrylate		bromo aromatic, high RI	03330-10	10 g	157
Triethylene glycol monomethyl ether monomethacrylate		hydrophilic	18556-500	500 g	157
3,3,5-Trimethylcyclohexyl methacrylate		aliphatic, bulky	02660-100	100 g	158
Undecyl methacrylate		hydrophobic	02544-25	25 g	158

## Adhesion Promoting Monomers

Functional groups known to increase adhesion of polymers to surfaces include phosphate and carboxylic acids (metal adhesion) and silyl ethers (glass/siliceous adhesion) which hydrolyze to give reactive Si-OH bonds. While these monomers are well studied examples, many monomers having functional groups such as acids, amines and hydroxyls can also impart polymer adhesion to various substrates. Please refer to other tables for a more expansive listing of these monomers. Also see Amine Monomer section.

	Polymerizable Sites	Polymerization Synthon	Add'l Reactive Functionality	Special Features			
Acrylic acid min. 99.5%	mono	acrylic	carboxylic acid	acid provides metal adhesion	00020-250	250 g	128
Bis(2-methacryloxyethyl) phosphate	dual	acrylic	phosphate	crosslinking monomer with adhesion promoting capabilities, good for metals	16041-10	10 g	132
4-Methacryloxyethyl trimellitic anhydride	mono	acrylic	anhydride	hydrolyzed acid offers improved adhesion	17285-10	10 g	148
3-Methacryloxypropyltrimethoxysilane	mono	acrylic	silyl ether	ethers react with siliceous surfaces to improve adhesion, glass pretreatment for polyacrylamide gels	02476-250	250 g	148
Monoacryloxyethyl phosphate	mono	acrylic	phosphate	used for introducing phosphorus into polymers, adhesion promoter	22468-10	10 g	150
Vinyltriethoxysilane	mono	vinyl	silyl ether	reactive silyl ethers hydrolyze affording bonding sites to siliceous surfaces	04537-50	50 g	160

## Amine Containing Monomers

Amines are among the most widely versatile functional groups. In biopolymers, amines are the key synthetic handle to build structure and architecture to a polymer. Amine groups can act as base catalysts, can be quaternized to yield aqueous soluble polymers and can function as ligands to a variety of metals. Amines are good nucleophiles and can be converted to a wide set of functional groups. Amines can form salts with carboxylic and phosphoric acids to form biologically interesting complexes and structures.

	Polymerizable Sites	Polymerization Synthon	Add'l. Reactive Functionality	Special Features	Catalog #	Size	Page
2-Acryloxyethyltrimethylammonium chloride	mono	acrylic	quaternary, HCl salt		17981-250	250 g	128
2-Aminoethyl methacrylate hydrochloride, min. 95%	mono	acrylic	primary, HCl salt		21002-10	10 g	129
N-(2-aminoethyl) methacrylamide hydrochloride	mono	acrylic	secondary, HCl salt		24833-5	5 g	129
N-(t-BOC-aminopropyl) methacrylamide	mono	acrylic	primary (protected)	solvent soluble, for primary amine polymers	24318-10	10 g	129
N-(3-Aminopropyl)methacrylamide hydrochloride, >98%	mono	acrylic	primary, HCl salt		21200-5	5 g	129
2-(tert-Butylamino)ethyl methacrylate, min. 90%	mono	acrylic	secondary	homopolymer T <sub>g</sub> = 33° C	01797-100	100 g	134
Diallylamine, min. 98%	dual	vinyl	secondary	forms cyclopolymers	21424-100	100 g	137
Diallyldimethylammonium chloride, 65% soln. in water	dual	vinyl	quat	cationic polymers via cyclopolymerization	15912-100	100 g	137
2-(N,N-Diethylamino)ethyl methacrylate	mono	acrylic	tertiary	homopolymer T <sub>g</sub> = 20° C	01872-500	500 g	137
2-(Diethylamino)ethylstyrene, mixed m,p-isomers, min. 95%	mono	vinyl	tertiary		24100-5	5 g	137
2-Diisopropylaminoethyl methacrylate	mono	acrylic	tertiary		24263-10	10 g	138
2-(N,N-Dimethylamino)ethyl acrylate	mono	acrylic	tertiary		02257-500	500 g	138
N-[2-(N,N-Dimethylamino)ethyl] methacrylamide	mono	acrylic	tertiary	hydrolytic stability	06172-5	5 g	138
N-[3-(N,N-Dimethylamino)propyl] acrylamide, min.95%	mono	acrylic	tertiary	homopolymer T <sub>g</sub> = 19° C	22018-10	10 g	138
N-[3-(N,N-Dimethylamino)propyl] methacrylamide	mono	acrylic	tertiary	hydrolytic stability	09656-100	100 g	139
3-Dimethylaminoneopentyl acrylate	mono	acrylic	tertiary		17970-10	10 g	139
2-(N,N-Dimethylamino)ethyl methacrylate, min. 99%	mono	acrylic	tertiary	homopolymer T <sub>g</sub> = 19° C	00213-500	500 g	138
Methacryloyl-L-Lysine	mono	acrylic	amino acid	zwitterionic, multifunctional, water soluble	24315-5	5 g	149
2-N-Morpholinoethyl acrylate, 95%	mono	acrylic	tertiary		17977-100	100 g	150
2-N-Morpholinoethyl methacrylate, 95%	mono	acrylic	tertiary		17978-100	100 g	150

## Crosslinking Acrylic Monomers - Difunctional

Difunctional monomers are useful for imparting crosslinking or branching sites to polymer architectures. The "spacer" group between the acrylic end groups often helps determine the physical and mechanical attributes of the resulting crosslinked polymer structure. Acrylic moieties are generally more reactive than methacrylic moieties and are thus used when faster reaction kinetics are desired (e.g. UV curable systems.)

	Add'l. Reactive Functionality	Special Features			
1-(Acryloyloxy)-3-(Methacryloyloxy)-2-Propanol	hydroxyl		25351-50	50 ml	129
Barium methacrylate, >95%	carboxylic acid	divalent metal atom salt	01994-50	50 g	130
2,2-Bis[4-(2-acryloxyethoxy)phenyl]propane	rigid hydrophobic, crosslinker	rigid, aromatic	04136-25	25 g	131
2,2-Bis[4-(2-hydroxy-3-methacryloxypropoxy)phenyl]propane		rigid, aromatic	03344-100	100 g	131
Bis(2-methacryloxyethyl) N,N'-1,9-nonylene bisacrylamate		flexible, hydrophobic	21619-50	50 g	131
Bis(2-methacryloxyethyl) phosphate	phosphate	adhesion promoter	16041-10	10 g	132
2,2-Bis(4-methacryloxyphenyl)propane		rigid, hydrophobic	01381-25	25 g	132
1,4-Butanediol diacrylate, min. 85%		aliphatic	02049-100	100 g	132
1,3-Butanediol dimethacrylate, 98%		aliphatic	02047-500	500 g	132
1,4-Butanediol dimethacrylate, min. 90%		aliphatic	05973-250	250 g	133
Copper (II) methacrylate	carboxylic acid	divalent metal atom salt	21222-25	25 g	135
trans-1,4-Cyclohexanediol dimethacrylate		rigid, aliphatic	18912-10	10 g	135
N,N'-Cystaminebisacrylamide, Electro Pure™	disulfide	reversible crosslinking	09809-5	5 g	135
1,10-Decanediol dimethacrylate		flexible, aliphatic, long chain	02140-25	25 g	136
1,4-Diacryloylpiperazine		rigid	21190-10	10 g	136
N,N-Diallylacrylamide	acrylic		01848-10	10 g	136
Diethylene glycol diacrylate		hydrophilic	02215-100	100 g	137
Diethylene glycol dimethacrylate		hydrophilic	02214-100	100 g	138
2,2-Dimethylpropanediol dimethacrylate		aliphatic	02276-100	100 g	139
N, N' Ethylene Bisacrylamide		hydrogel crosslinker	09811-1	1 g	140
Ethylene glycol diacrylate		UV curable, acrylate	02302-25	25 g	140
Ethylene glycol dimethacrylate, 99.7%		hydrophilic, high purity	24896-250	250 g	141
Ethylene glycol dimethacrylate, min 88%		hydrophilic	02303-250	250 g	141
Ethylene glycol dimethacrylate, min 98%		hydrophilic	24030-250	250 g	141
Fluorescein dimethacrylate		fluorescent	23589-100	100 mg	141
N,N'-Hexamethylenebisacrylamide		hydrolytic stability	01495-5	5 g	143
1,6-Hexanediol diacrylate		aliphatic	23671-100	100 g	144
1,6-Hexanediol dimethacrylate, min 98%		aliphatic	23672-100	100 g	144
2,2-Bis[4-(2-hydroxy-3-methacryloxypropoxy)phenyl]propane		rigid, aromatic	03344-500	500 g	144
Lead acrylate			02968-100	100 g	146
Magnesium acrylate	carboxylic acid	divalent metal atom salt	02467-10	10 g	147
N,N'-Methylenebisacrylamide, Chemzymes®, Ultra Pure, Purity >99%		hydrogel crosslinker	00719-100	100 g	149
1,9-Nonanediol dimethacrylate		flexible, aliphatic	00801-10	10 g	151
1,5-Pentanediol dimethacrylate		aliphatic	04260-25	25 g	152
1,4-Phenylene diacrylate		rigid, aromatic	06389-10	10 g	153
Tetraethylene glycol dimethacrylate		aliphatic	02654-50	50 g	156
Tricyclodecane dimethanol diacrylate			25110-50	50 g	157
Triethylene glycol diacrylate		UV curable, acrylate > methacrylate	02655-250	250 g	157
Triethylene glycol dimethacrylate			01319-250	250 g	157
Triethylene glycol dimethacrylate, min 95%		aliphatic	24034-100	100 g	157
Zinc (di)methacrylate	carboxylic acid	ionomeric crosslinker	03011-100	100 g	160

## Crosslinking Acrylic Monomers - Multifunctional

Typically used for generating highly crosslinked polymer structures, these monomers increase polymer toughness, modulus and solvent resistance. For UV curable formulations, multifunctional acrylates are typically faster reacting than their methacrylate analogs.

	Polymerizable Sites	Add'l. Reactive Functionality	Special Features			
Dipentaerythritol pentaacrylate (mixture of tetra-, penta-, hexaacrylate)	multiple		high crosslinking efficiency	16311-500	500 g	139
Pentaerythritol tetraacrylate	multiple			01547-100	100 g	152
Pentaerythritol triacrylate	tri	hydroxyl	multifunctional crosslinker	04259-100	100 g	152
1,1,1-Trimethylolpropane triacrylate	tri		useful for UV cure	02658-250	250 g	158
1,1,1-Trimethylolpropane trimethacrylate	tri		useful for UV cure	02659-250	250 g	158

## Dual Reactive Acrylic Monomers

It is often desirable to synthesize polymer architectures that are capable of further reaction to incorporate new functionality, graft new polymer chains, attach drugs or biomolecules, or make the polymer respond intelligently to changes in its environment. This diverse set of monomers have easily polymerizable carbon-carbon double bonds yet contain a secondary reactive group that can be elaborated in a multitude of ways. Some reactive groups e.g., carboxylic acid in o-nitrobenzyl methacrylate (*Cat. #24360*) are masked and are revealed by simple deprotection schemes.

	Add'l. Reactive Functionality	Special Features			
Acrylic anhydride min. 90%	anhydride (protected)	allows formation of cyclic anhydrides	00488-50	50 g	128
4-(2-Acryloxyethoxy)-2-hydroxybenzophenone	phenol	benzophenone is a strong UV absorber, potential UV radical initiator	19931-10	10 g	128
N-Acryloxysuccinimide	ester (protected)	derivitization through ester link using mild reaction conditions, good for biologically active compounds	19930-1	1 g	128
Allyl methacrylate	olefin	allyl double bond is less reactive than methacrylate, good for post polymerization reactions	01643-500	500 g	129
2-Aminoethyl methacrylate hydrochloride, min. 95%	amine (as HCl salt)	can be used to prepare polymers with amine functionality	21002-10	10 g	129
N-(t-BOC-aminopropyl)methacrylamide	amine (protected)	t-BOC easily deprotected to yield primary amine	24318-10	10 g	129
N-(3-Aminopropyl)methacrylamide hydrochloride, >98%	amine (as HCl salt)	can be used to prepare polymers with amine functionality	21200-5	5 g	129
2-Bromoethyl acrylate, min. 95%	halogen	allows synthesis of heavy atom polymers, Br can be displaced with various nucleophiles	02015-10	10 g	132
N-(iso-Butoxymethyl)methacrylamide, min. 95%	oxo-methylene	acid + heat catalyzed crosslinking monomer	19221-100	100 g	133
tert-Butyl acrylate	carboxylic acid (protected)	used in photoresist formulations	02039-250	250 g	133
2-(tert-Butylamino)ethyl methacrylate, min. 90%	amine (secondary)	builds polymers with secondary amines, potential H-Bonding sites	01797-100	100 g	134
2-Carboxyethyl acrylate, 90%	carboxylic acid	lower volatility analog to acrylic acid	16719-100	100 g	134
2-Cinnamoyloxyethyl acrylate	$\alpha$ , $\beta$ unsat. ester	photocrosslinking monomer	24014-10	10 g	135
Cinnamyl methacrylate	functionalized styrene	photocrosslinking monomer	02092-5	5 g	135
2-Cyanoethyl acrylate	nitrile	polar building block monomer, can use nitrile for functionalization	01829-100	100 g	135
N,N-Diallylacrylamide	allyl	crosslinking monomer	01848-10	10 g	136
Dicyclopentenylacrylate	ethenyl	endo cyclic olefin does not readily polymerize, can be post reacted e.g. oxidative crosslinking	15797-25	25 g	137
Glycerol monomethacrylate, mixture of isomers	hydroxyl	isomer mixture, useful in hydrogel preparation, -OH can be functionalized in multiple ways	04180-25	25 g	141

			Catalog #	Size	Page
Glycidyl methacrylate	Add'l. Reactive Functionality epoxide	Special Features versatile, introduces reactive sites into polymers, can be derivatized in post polymerization reactions with various nucleophiles	02607-500	500 g	142
2-Hydroxyethyl acrylate	hydroxyl	hydrophilic monomer building block	01902-250	250 g	144
Hydroxypropyl methacrylate, mixture of isomers	hydroxyl	hydrophilic monomer building block, can be post reacted through hydroxyl	00730-500	500 g	146
N-(2-Hydroxypropyl)methacrylamide	hydroxyl	hydrophilic monomer building block, can be post reacted through hydroxyl group	08242-10	10 g	145
3-Methacryloxypropyltrimethoxysilane	silyl ether	ethers react with glass surfaces to improve adhesion, glass pretreatment for polyacrylamide gels	02476-250	250 g	148
Methacryloyl chloride, min. 80%	halogen	reactive building block monomer. Also used to prepare other specialty monomers	01518-50	50 g	148
Methacryloyl fluoride	halogen	reactive building block monomer. Also used to prepare other specialty monomers	17414-25	25 g	149
Methyl 2-cyanoacrylate	cyano	used in rapidly setting cements and adhesives	01520-10	10 g	149
N-Methylolacrylamide, 48% soln. in water		acid catalyzed post polymerization crosslinking capabilities	02518-250	250 g	149
Monoacryloxyethyl phosphate	phosphate	used for introducing phosphorus into polymers, adhesion promoter	22468-10	10 g	150
o-Nitrobenzyl methacrylate, min. 95%	carboxylic acid (protected)	nitrobenzyl ester removed by UV irradiation to afford carboxylic acid, used for catalysis, photoresists, latent reactive acid	24360-10	10 g	151
N-(Phthalimidomethyl)acrylamide	active methylene	used in photoresist formulations	19390-25	25 g	153
Propargyl acrylate	acetylenic	acetylenic group can be post functionalized, possible oxidative crosslinking monomer	02964-25	25 g	154
Propargyl methacrylate	acetylenic	acetylenic group can be post functionalized, possible oxidative crosslinking monomer	02965-25	25 g	154
Tricyclodecane dimethanol diacrylate			25110-50	50 g	157

## Epoxides / Anhydrides / Imides

Unlike ethylenic monomers which polymerize through free radical processes, epoxide resins react with amines, carboxylic acids, anhydrides, etc. to form polymers displaying a range of characteristics from tough and durable to soft and adhesive.

	Polymerizable Sites	Special Features			
Allyl glycidyl ether	mono	can react either epoxide or allyl group into polymers	19191-50	50 g	129
4,4'-Bisphenol A Bis-(N-Methylphthalimide)	dual	useful monomer for synthesis of high temp. polymers	24284-25	25 g	132
4,4'-Bisphenol A Dianhydride	dual	useful monomer for synthesis of high temp. polymers	24283-25	25 g	132
1,4-Butanediol diglycidyl ether	dual	can be used for post polymerization crosslinking of amine polymers	01795-50	50 g	132
Citraconic anhydride	mono	reversible protection of amino groups in protein synthesis	01824-100	100 g	135
Glycidyl butyl ether	mono	used in epoxy reactions as monofunctional diluent	05678-500	500 g	141
Glycidyl cinnamate	mono	polymerizable with various carboxylic acids, alcohols and amines; Olefin moiety can be photocrosslinked	16090-10	10 g	142
Glycidyl Glycerol-Ether, Polyfunctional	multiple	isomer mixture, efficient crosslinking monomer for various nucleophilic agents	09221-50	50 g	142

# Monomers

			Catalog #	Size	Page
Glycidyl nonylphenyl ether	Polymerizable Sites mono	Special Features epoxy functional	18236-100	100 g	142
Methacrylic anhydride		uncrosslinked	01517-100	100 ml	147
Ethylene glycol diglycidyl ether (Quetol 651)	dual	active crosslinking monomer for active nucleophiles, hydrophilic, used for hydrogels	01479-100	100 g	141
Propylene Oxide, EM Grade	mono	basic bulding block monomer for water dispersible polymers	00236-1	1 pint	155
Propylene glycol diglycidyl ether	dual	versatile crosslinker for amine, hydroxyl and carboxylate systems	24044-100	100 g	155
Triglycidyl isocyanurate	tri	crosslinking monomer for epoxy, urethane systems	16173-50	50 g	158

## Fluorescent Acrylic Monomers

Monomers with fluorescent tags are often used to build polymers that can be detected at very low concentrations using fluorescence spectroscopy. Polymer migration and diffusion has been studied using fluorescent tags. Polymer microspheres containing fluorescent groups are used routinely for flow cytometry and medical diagnostic assays.

	Monomer Type	Special Features			
9-Anthracenylmethyl methacrylate	neutral	Ex. max = 362 nm, Em. max = 407 nm, yellow crystals	23587-100	100 mg	130
	neutral	Ex. max = 362 nm, Em. max = 407 nm, yellow crystals	23587-1	1 g	130
3,8-dimethacryloyl ethidium bromide	ionic	Ex. max = 439 nm, Em. max = 512 nm, insoluble in water	23590-100	100 mg	138
Fluorescein dimethacrylate	neutral	Ex. max = 470 nm, Em. = 511 nm	23589-100	100 mg	141
Methacryloxyethyl thiocarbamoyl rhodamine B	ionic, mixture isomers	Ex. max = 548 nm, Em. max = 570 nm, purple crystals	23591-100	100 mg	147
2-Naphthyl acrylate	neutral	hydrophobic, fluorescent monomer	06024-1	1 g	150
2-Naphthyl methacrylate	neutral	Ex. max = 285 nm, Em. min = 345 nm	23602-100	100 mg	150
Nile Blue Acrylamide		dark blue crystals; Em. Max. 674 nm in H <sub>2</sub> O	25395-100	100 mg	150
Nile Blue Methacrylamide		dark purple crystals; Em. Max. 674 nm in H <sub>2</sub> O	25355-100	100 mg	151
O-Methacryloyl Hoechst 33258	neutral	Ex. max = 355 nm, Em. max = 497 nm, off-white crystals	23592-100	100 mg	148
1-Pyrenylmethyl methacrylate	neutral	Ex. max = 339 nm, Em. min = 394 nm, pale yellow crystals	23588-100	100 mg	155

## Fluorinated Acrylic Monomers

Monomers containing fluorine provide polymers with unique low energy surfaces. Materials made from these monomers are typically chemical resistant and very hydrophobic.

	Homopolymer Tg (°C)			
Bis-(2,2,2-Trifluoroethyl) Itaconate		24972-25	25 g	131
Bis-(1,1,1,3,3,3-Hexafluoroisopropyl) Itaconate		24971-25	25 g	131
1H,1H,7H-Dodecafluoroheptyl methacrylate	13	00767-25	25 g	140
1H,1H,2H,2H-Heptadecafluorodecyl acrylate	0	19227-100	100 g	142
1H,1H,2H,2H-Heptadecafluorodecyl methacrylate	40	19226-25	25 g	143
1H,1H-Heptafluorobutyl acrylate	-30	21039-25	25 g	143
Hexafluoro-iso-propyl methacrylate		02401-10	10 g	143
1H,1H,3H-Hexafluorobutyl acrylate	-22	05631-10	10 g	143
1H,1H,3H-Hexafluorobutyl methacrylate		05632-10	10 g	143
1,1,1,3,3,3-Hexafluoroisopropyl Acrylate		24970-25	25 g	143
1H,1H,5H-Octafluoropentyl acrylate	-35	21044-25	25 g	151
1H,1H,5H-Octafluoropentyl methacrylate, min. 98%	36	21045-25	25 g	151
Pentafluorophenyl acrylate		06349-5	5 g	152
Pentafluorophenyl methacrylate, 95%		06350-5	5 g	152
1H,1H,3H-Tetrafluoropropyl methacrylate		07577-25	25 g	156
2,2,2-Trifluoroethyl acrylate	-10	01718-25	25 g	157
2,2,2-Trifluoroethyl methacrylate	80	02622-25	25 g	157

## High/Low Refractive Index Monomers

Polymers that interact with light can be modified by optimizing their refractive index properties. These materials are useful in many types of optical applications: lenses, optical switches, optical fiber coatings, etc. Low RI monomers are typically highly fluorinated. Polymers using perfluorinated monomers are often hydrophobic and exhibit very low surface energies.

	Polymerizable Sites	Homopolymer Tg (°C)	Special Features			
Allyl Phenyl Ether, 98%	vinyl		hydrophobic, high RI (ca 1.52)	24894-100	100 g	129
Benzyl acrylate, ~99%	acrylic	6	high RI (ca 1.55), non-halogenated	01997-100	100 g	130
Benzyl methacrylate, min. 95%	acrylic	54	high RI (ca 1.57), non-halogenated	02000-100	100 g	130
N-Benzylmethacrylamide	acrylic		high RI (ca 1.60), non-halogenated	17969-25	25 g	130
2-(9H-Carbazole-9-yl) ethyl methacrylate	acrylic		high RI (ca 1.69)	24372-1	1 g	134
4-Chlorophenyl acrylate	acrylic	58	high RI (ca 1.55), chlorinated aromatic	01331-10	10 g	134
3-Chlorostyrene, 98%		90	high RI (ca 1.57)	02042-5	5 g	135
1H,1H,7H-Dodecafluoroheptyl methacrylate	acrylic	13	low RI (ca 1.36), fluorinated aliphatic	00767-25	25 g	140
1H,1H,2H,2H-Heptadecafluorodecyl acrylate	acrylic	0	low RI (ca 1.34), fluorinated	19227-25	25 g	142
1H,1H,2H,2H-Heptadecafluorodecyl methacrylate	acrylic	40	low RI (ca 1.35), fluorinated	19226-25	25 g	143
1H,1H-Heptafluorobutyl acrylate	acrylic	-30	low RI (ca 1.37), fluorinated	21039-25	25 g	143
Hexafluoro-iso-propyl methacrylate	acrylic		low RI (ca 1.38), fluorinated	02401-10	10 g	143
1H,1H,3H-Hexafluorobutyl acrylate	acrylic	-22	low RI (ca 1.39), fluorinated	05631-10	10 g	143
1H,1H,3H-Hexafluorobutyl methacrylate	acrylic		low RI (ca 1.40), fluorinated	05632-10	10 g	143
1H,1H,5H-Octafluoropentyl acrylate	acrylic	-35	low RI (ca 1.38), fluorinated	21044-25	25 g	151
1H,1H,5H-Octafluoropentyl methacrylate, min. 98%	acrylic	36	low RI (ca 1.39), fluorinated	21045-25	25 g	151
Pentabromophenyl acrylate	acrylic		high RI (ca 1.7), brominated aromatic	06344-10	10 g	152
Pentabromophenyl methacrylate	acrylic		high RI (ca 1.7), brominated aromatic	04253-5	5 g	152

# Monomers

				Catalog #	Size	Page
1H,1H,3H-Tetrafluoropropyl methacrylate	Polymerizable Sites: acrylic	Homopolymer Tg (°C):	Special Features: low RI (ca 1.4), fluorinated aliphatic	07577-25	25 g	156
2,4,6-Tribromophenyl acrylate	acrylic		high RI (ca 1.6), brominated aromatic	03330-10	10 g	157
2,2,2-Trifluoroethyl acrylate	acrylic	-10	low RI (ca 1.44), fluorinated aliphatic	01718-25	25 g	157
2,2,2-Trifluoroethyl methacrylate	acrylic	80	low RI (ca 1.41), fluorinated aliphatic	02622-25	25 g	157
N-Vinylcarbazole	vinyl		high RI (ca 1.68)	02429-25	25 g	159

## Hydroxy Containing Monomers

Hydroxyl groups have utility as hydrogen bonding sites and can provide polymers with compatibility for water or polar solvents. These versatile functional groups can be derivitized broadly. Polymers containing free –OH groups can be post reacted with acids, epoxies, isocyanates, etc. to create novel polymer properties and architectures.

	Polymerizable Sites	Polymerizable Synthon	Add'l. Reactive Functionality	Special Features			
4-(2-Acryloxyethoxy)-2-hydroxybenzophenone	mono	acrylic	phenol	UV absorbing, aromatic	19931-10	10 g	128
4-tert-Butoxystyrene, min. 99%	mono	styrene	phenol (protected)	deprotect with heat	21760-10	10 g	133
Glycerol monomethacrylate, mixture of isomers	mono	acrylic	hydroxyl	useful for hydrogels	04180-25	25 g	141
(HEMA 10) Poly Ethoxy (10) ethyl methacrylate	mono	acrylic	hydroxyl	hydrophilic	24890-100	100 g	142, 153
2-Hydroxyethyl acrylate	mono	acrylic	hydroxyl	homopolymer Tg = 15° C	01902-250	250 g	144
Hydroxypolyethoxy (10) Allyl Ether, 98%	mono	vinyl	hydroxyl	hydrophilic	24899-100	100 g	145
N-(2-Hydroxypropyl) methacrylamide	mono	acrylic	hydroxyl	hydrolytic stability	08242-10	10 g	145
Hydroxypropyl methacrylate, mixture of isomers	mono	acrylic	hydroxyl	aliphatic	00730-1	1 kg	146
Poly(propylene glycol) (300) monomethacrylate	mono	acrylic	hydroxyl	slightly hydrophilic, long chain graft	15934-250	250 g	154
Sodium 1-Allyloxy-2-hydroxypropyl Sulfonate				sulfonate salt	24898-100	100 g	155
1,1,1-Trimethylolpropane diallyl ether	mono	vinyl	hydroxyl	hydrophilic crosslinker	05500-50	50 g	158
1,1,1-Trimethylolpropane monoallyl ether	mono	vinyl	hydroxyl	hydrophilic crosslinker	15914-50	50 g	158

## Mono and Difunctional Glycol Oligomeric Monomers

Ethylene glycol units are strongly hydrophilic through their multiple H-bonding sites. Monomers of this type are useful in the construction of hydrogels and water compatible polymer structures. New research suggests that bioactive molecules e.g. drugs with attached PEG chains have improved bioavailability characteristics.

	Polymerizable Synthon (Sites)	Add'l Reactive Group	Approx. # of Glycol Groups	Approx. MW.	Special Features	Catalog #	Size	Page
Poly(ethylene glycol) diacrylate	acrylic (2)	none	4 to 6	314	hydrophilic crosslinking monomer	00669-250	250 g	153
	acrylic (2)	none	8 to 10	500		01871-250	250 g	153
	acrylic (2)	none	90 to 100	4,100		15246-1	1 g	153
Poly(ethylene glycol) dimethacrylate	methacrylic (2)	none	4 to 6	350	hydrophilic crosslinking monomer	00096-100	100 g	154
	methacrylic (2)	none	8 to 10	550		15179-100	100 g	154
	methacrylic (2)	none	13 to 15	750		02364-100	100 g	154
	methacrylic (2)	none	21 to 25	1,150		15178-100	100 g	154
Poly(ethylene glycol) monomethacrylate	methacrylic (1)	hydroxy	4 to 6	270	long chain monomer, more polar than PO analog, can post react -OH to add additional functionality	16712-100	100 g	154
	methacrylic (1)	hydroxy	8 to 10	480		16713-100	100 g	154
Poly(ethylene glycol) monomethyl ether monomethacrylate	methacrylic (1)	none	4 to 6	280	adds hydrophilic grafts to polymers; adds long chain hydrophilic graft to polymer chain	16664-100	100 g	124
	methacrylic (1)	none	8 to 10	490		16665-100	100 g	124
	methacrylic (1)	none	21 to 25	1,090		16666-100	100 g	124
Poly(propylene glycol) monomethacrylate	methacrylic (1)	hydroxy	3 to 5	440	propylene glycol is more hydrophobic than ethylene glycol long chain monomer, can post react -OH to add functionality	15934-250	250 g	154
Poly(propylene glycol) dimethacrylate	methacrylic (2)	none	5 to 7	550	less polar than EO analog	04380-250	250 g	154

## Styrenic Monomers

Popular alternatives to acrylic and related monomers, styrenic monomers generally provide polymers of higher glass transition temperature, higher modulus, increased hydrophobic character and nominal UV absorbance. As such, coatings made with high concentrations of styrenic monomers can yellow with time if exposed to UV light. Crosslinked styrene resins (especially in microsphere form) are tough and chemically resistant. These form the basis for ion exchange resins and microbeads used as supports for biochemical reactions.

	Monomer Type	Add'l. Reactive Functionality	Homopolymer Tg (°C)	Special Features			
4-tert-Butoxystyrene, min. 99%	protected group	phenol		deprotection gives phenol	21760-10	10 g	133
4-(tert-Butyl)styrene, 95%	neutral		132	high Tg monomer	02606-25	25 g	134
4-Chloromethylstyrene	reactive	halogen		Merrifield resin building block	22193-25	25 g	134
Chloromethylstyrene, 43% para, 57% meta, ~96%	reactive	halogen		Merrifield resin building block	02718-500	500 g	134
3-Chlorostyrene, 98%			90		02042-5	5 g	135
2-(Diethylamino)ethylstyrene, mixed m,p-isomers, min. 95%					24100-5	5 g	137
Divinylbenzene (50 - 60% Active)					01892-100	100 g	140
Divinylbenzene (80% Active)	neutral	vinyl		aromatic, rigid crosslinker	22478-100	100 g	140
4-Methylstyrene	neutral		108	aromatic, hydrophobic	04234-100	100 g	150
2-Methylstyrene	neutral		120	aromatic, hydrophobic	04581-5	5 g	149
4-Nitrostyrene	reactive	nitro		nitro can be reduced to amine	02634-5	5 g	151
Styrene	neutral		100		00660-500	500 g	156
4-Vinylbenzoic acid	reactive	carboxylic acid		versatile -COOH synthesis handle	04485-5	5 g	159

## Sulfonate Containing Monomers

	Polymerizable Sites	Polymerization Synthon	Form of Acid	Special Features			
2-Sulfoethyl methacrylate, >90%	mono	acrylic	sulfonic acid	water soluble	02597-50	50 g	156
3-Sulfopropyl methacrylate, potassium salt, 98%	mono	acrylic	sulfonic acid salt	water soluble	17210-100	100 g	156
3-Sulfopropyl dimethyl-3-methacrylamidopropyl ammonium, inner salt				Zwitterionic detergent and crosslinker.	16570-5	5 g	156

## UV (light) Active Monomers

Polymers with aliphatic backbones often show little absorbance of light and usually do not absorb in the near and mid UV spectral range. UV absorbing monomers improve the capture of light at these wavelengths. These absorbers can be used to shield the polymer system or an underlying substrate from degradation by UV light, e.g. phenethyl methacrylate containing polymers for optical lenses. Additionally, some UV absorbing materials can act as sensitizers to promote photochemical reactions.

	Polymerization Synthon	Add'l. Reactive Functionality	Homopolymer Tg (°C)	Special Features			
2(5-Benzotriazolyl 3-t-butyl 4-methacryloxy)phenyl propionic acid, isooctyl ester	acrylic			UV absorber	24241-1	1 g	130
2-Cinnamoyloxyethyl acrylate	acrylic	$\alpha,\beta$ unsat.		photo reactive ester	24014-10	10 g	135
Cinnamyl methacrylate	acrylic	$\alpha,\beta$ unsat.		photo reactive ester	02092-5	5 g	135
Glycidyl cinnamate	epoxide	ethenyl		photo reactive	16090-10	10 g	142

					Catalog #	Size	Page
	Polymerization Synthon	Add'l. Reactive Functionality	Homopolymer Tg (°C)	Special Features			
2-(2'-Methacryloxy-5'-methylphenyl) benzotriazole	acrylic			UV absorber	21871-25	25 g	148
2-(2'-Hydroxy-5'-methylphenyl)-benzotriazole	none			UV absorber	21363-25	25 g	144
Phenyl acrylate, min. 95%	acrylic		57	UV absorbing	02642-10	10 g	153
Phenyl methacrylate, >95%	acrylic		110	moderate UV absorbing	02644-10	10 g	153
2-Phenylethyl acrylate, min. 92%	acrylic		-3	moderate UV absorbing, aromatic, hydrophobic	02834-100	100 g	153
2-Phenylethyl methacrylate, min. 92%	acrylic		26	moderate UV absorbing	02911-100	100 g	153
N-(Phthalimidomethyl)acrylamide	acrylic			used in photoresist formulations	19390-25	25 g	153

## Vinyl and Ethenyl Monomers




An alternative choice to styrenics and acrylics, these monomers are often used to create polymers with inert main chain features. Many of these monomers can be polymerized via metallocene or other metal mediated polymerization processes.

	Polymerizable Sites	Add'l. Reactive Functionality	Special Features			
Acrolein, distilled, 99% EM grade	mono	aldehyde		00016-5	5x10 ml	128
Allyl glycidyl ether	mono	epoxide		19191-50	50 g	129
2,2-Bis[4-(2-acryloxyethoxy)phenyl] propane	dual	rigid crosslinker		04136-25	25 g	131
3-Butene-1,2,3-tricarboxylic acid	mono			01396-25	25 g	133
Diallyl Maleate	dual vinyl	vinyl	hydrophilic	02156-250	250 g	136
Diallyl Maleate - 99% Active	dual vinyl	vinyl	hydrophilic	24892-100	100 g	136
Diallyl phthalate, practical	dual			02159-1	1 kg	136
Diallylamine, min. 98%	dual	amine (secondary)		21424-100	100 g	137
Diallyldimethylammonium chloride, 65% soln. in water	dual			15912-100	100 g	137
$\alpha,\alpha$ -Dimethyl-3-isopropenylbenzyl isocyanate	mono	isocyanate		19706-100	100 g	138
Divinyl Adipate	dual			25338-100	100 g	139
Divinyl glycol	dual			18252-25	25 g	139
Divinyl sebacate	dual			04632-5	5 g	139
N-Methyl N-vinyl acetamide	mono		hydrophilic	22065-25	25 g	149
2-Propene-1-sulfonic acid, sodium salt, 35% soln. in water	mono		sulfonic acid, Na salt	00064-10	10 g	154
Triallyl cyanurate	tri			01236-100	100 g	157
Triethylene glycol divinyl ether	dual			19560-100	100 g	157
1,1,1-Trimethylolpropane diallyl ether	multi	hydroxyl		05500-50	50 g	158
1,1,1-Trimethylolpropane monoallyl ether	mono	hydroxyl		15914-50	50 g	158
N-Vinyl acetamide (NVA)	mono			24806-50	50 g	158
Vinyl azlactone	mono			21329-5	5 g	159
Vinyl benzoate	mono		aromatic	02664-10	10 g	159
Vinyl butyrate	mono			02835-10	10 g	159
Vinyl 2-furoate	mono			02829-10	10 g	159
Vinyl methyl sulfone	mono			17206-5	5 g	126

				Catalog #	Size	Page
Vinyl octadecyl ether	Polymerizable Sites: mono	Add'l. Reactive Functionality:	Special Features: hydrophobic	01728-100	100 g	159
Vinyl stearate, min. 95%	mono			01784-100	100 g	159
N-Vinyl-2-pyrrolidone	mono		hydrophilic	04000-250	250 g	160
N-Vinylcaprolactam, min. 98%	mono			16818-10	10 g	159
N-Vinylcarbazole	vinyl		high RI (ca 1.68)	02429-25	25 g	159
Vinylferrocene	mono		iron complex	04503-1	1 g	159
1-Vinylimidazole	mono			01726-100	100 g	160
4-Vinylpyridine	mono			02668-100	100 g	160
Vinyltriethoxysilane	mono	silyl ether	adhesion promoter	04537-50	50 g	160
Vinyltriphenylsilane				21543-10	10 g	160

## Polymerization Inhibitors for Monomers

In the Monomer section, we list monomers which are inhibited with a variety of polymerization inhibitors. These inhibitors are chosen for effectiveness, and minimum color formation on storage. Below is a table identifying the inhibitors used and their structures.

Inhibitor	Chemical Name	Structure	Inhibitor	Chemical Name	Structure
HQ	Hydroquinone		Triethylene diamine (DABCO®)	2,4-Diazabicyclo[2.2.2]octane	
MEHQ	Hydroquinone monomethyl ether		t-Butylcatechol	4-t-butylcatechol	
PTZ	Phenothiazine		BHT(butylated hydroxytoluene)	2,6-di-t-Butyl-4-methyl-phenol	
Methylene Blue	3,7-Bis(dimethylamino)-phenazathionium chloride		p-t-Butylphenol	4-t-Butylphenol	

### Corporate Headquarters

**Polysciences, Inc.**  
400 Valley Road  
Warrington, PA 18976

(800) 523-2575 / (215) 343-6484  
(800) 343-3291 / (215) 343-0214 fax  
info@polysciences.com

### European Office

**Polysciences Europe GmbH**  
Handelsstrasse 3  
D-69214 Eppelheim, Germany

(49) 6221-765767  
(49) 6221-764620 fax  
info@polysciences.de

### Asia Pacific Office

**Polysciences Asia Pacific, Inc.**  
2F-1, 207 Tunhwa North Road  
Taipei, Taiwan 10595

(886) 2 8712 0600  
(886) 2 8712 2677 fax  
info@polysciences.tw