

## TECHNICAL DATA SHEET 558

Page 1 of 1

# BioMag® Goat anti-Mouse IgM

## Description

BioMag Goat anti-Mouse IgM is a suspension of BioMag particles approximately 1.6 $\mu$ m in size, which are covalently coated with goat anti-Mouse IgM antibody. The suspension is supplied in phosphate buffered saline (pH 7.4) with EDTA and sodium azide added. After shaking vigorously or vortexing, BioMag Goat anti-Mouse IgM is ready to use.

## Characteristics

Mean Diameter: ~1.6 $\mu$ m  
Particle Concentration: 1 mg/ml  
Binding Capacity: 1mL (1mg) of BioMag Goat anti-Mouse IgM will bind > 0.2mg of mouse IgM.

## Procedure

Researchers are advised to optimize the use of BioMag in any application as procedures designed by other manufacturers may not be ideal.

BioMag Goat anti-Mouse IgM is suitable for use as a secondary antibody in radioimmunoassays and enzyme immunoassays that utilize a mouse IgM primary monoclonal antibody. In cell sorting methods utilizing a mouse monoclonal (IgM) antibody against a specific cell population, BioMag Goat anti-Mouse IgM may be used to separate the cells of interest from a cell mixture. (For additional information on the use of this product in cell separation, please see Technical Data Sheet 528, 'BioMag and Cell Sorting.')

Separation may be accomplished magnetically.

## Storage and Stability

Store at 4°C. Freezing, drying, or centrifuging BioMag may result in irreversible aggregation and loss of binding activity.

## Safety

This particle suspension contains sodium azide. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Upon disposal of material, flush with a large volume of water to prevent azide accumulation. Please consult the Material Safety Data Sheet for more information.

**This product is for research use only and is not intended for use in humans or for *in vitro* diagnostic use.**

## Ordering Information

| Cat. #    | Description               | Size  |
|-----------|---------------------------|-------|
| 84350-50  | BioMag® Goat anti-Rat IgM | 50ml  |
| 84350-500 | BioMag® Goat anti-Rat IgM | 500ml |

### To Order

In The U.S. Call: 1-800-523-2575 • 215-343-6484  
In The U.S. FAX: 1-800-343-3291 • 215-343-0214

In Germany Call: (49) 6221-765767  
In Germany FAX: (49) 6221-764620

Order online anytime at [www.polysciences.com](http://www.polysciences.com).

Should any of our materials fail to perform to our specifications, we will be pleased to provide replacements or return the purchase price. We solicit your inquiries concerning all needs for life sciences work. The information given in this bulletin is to the best of our knowledge accurate, but no warranty is expressed or implied. It is the user's responsibility to determine the suitability for their own use of the products described herein, and since conditions of use are beyond our control, we disclaim all liability with respect to the use of any material supplied by us. Nothing contained herein shall be construed as a recommendation to use any product or to practice any process in violation of any law or any government regulation.