



U.S. Corporate Headquarters  
400 Valley Road  
Warrington, PA 18976  
1(800) 523-2575 / (215) 343-6484  
Fax: 1(800)343-3291  
Email: info@polysciences.com  
www.polysciences.com

Polysciences Europe GmbH  
Handelsstrasse 3 D-69214  
Eppelheim, Germany  
+(49) 6221-765767  
Fax: +(49) 6221-764620  
Email: info@polysciences.de

## TECHNICAL DATA SHEET 869

Page 1 of 1

# PC800 A&B

## Versatile, Fast Curing Acrylic Adhesive

### DESCRIPTION

**PC800 A&B** is a fast curing versatile acrylic adhesive. The product has excellent adhesion to many metals, plastics and composites. It is designed for applications needing good fast setting bonds at room temperature

### UNCURED PROPERTIES

	Part A	Part B	Mix
Type	Acrylic	Acrylic	Acrylic
Color	Colorless	Yellow	Colorless
Viscosity	Gel	Gel	Gel

### PROCESS PARAMETERS

The product is best dispensed through a cartridge static tip mixer.

#### Mix Ratio:

	Part A	Part B
By Weight	1	1
By Volume	1	1

	Mixed Product
Pot Life @25°C	10 minutes
Bond Set Time @25°C	30 minutes
Cure Time @25°C	8 hours

### CURED PROPERTIES

Property	Units	Value
Hardness @25°C	Shore D	65
Hardness @100°C	Shore D	40
Adhesion to Aluminum	psi	2200
Peel Strength	pli	25
Temperature Range of Use	°C	-50 to 120

### STORAGE AND HANDLING

#### Shipment

Room temperature

#### Storage

Room temperature

#### Safety

Refer to MSDS for details

All values are considered typical based on tests believed to be accurate. Polysciences, Inc. may change the data as appropriate.

### TO ORDER

**In The U.S. Call:** 1(800) 523-2575 • (215) 343-6484  
**In The U.S. Fax:** 1(800) 343-3291 • (215) 343-0214  
**In Germany Call:** +(49) 6221-765767  
**In Germany Fax:** +(49) 6221-764620

Order online anytime at [www.polysciences.com](http://www.polysciences.com)