

# Contract Packaging & Filling Capabilities

*Polysciences Inc. has conducted filling and private labeling for large Fortune 500 clients since the 1990s. Contact us to discuss your unique packaging requirements.*



## Clean Room Filling

- Explosion proof and Class 100 rated clean room
- Fill quantities from a few grams to 5 gallon pails
- Ultra Pure (USP grade or better) water solutions & sterile alcohol filling



## Specialized Small Packaging

- Inert gas (nitrogen or argon) filling
- Class 100,000 rated c-GMP areas equipped with glove boxes
- Manual sealing of glass ampoules containing air & moisture sensitive chemicals



## Filling Stations Large Quantity Packaging

- Packette filling of solid or liquid, ~70mm to ~14cm wide
- Non-flammable cold and hot filling of tubes, 10ml to 8oz
- Single & multi-head solid or liquid filling with one or two sided labeling, up to 1 gallon containers
- Filling capability up to tote quantities



For more information please visit us at [www.polysciences.com/confill](http://www.polysciences.com/confill)