

TECHNICAL DATA SHEET 404

Page 1 of 1

Polyballs

Description

Polyballs range in size from 1/16" (~1500µm) to 3/8" (~9.5mm). Some polystyrene-based Polyballs are available in functionalized or dyed versions, or with special surfaces for biological applications.

Sizes smaller than 100µm are available as Polybeads®. Please ask for Technical Data Sheet #238 for more details. For particles between 105-500µm or for glass beads, please see our Microspheres Catalog or check our website at www.polysciences.com.

Storage and Stability

Store dry product at room temperature.

This product is for research use only and is not intended for use in humans or for *in vitro* diagnostic use.

Ordering Information

Cat. #	Description	Diameter	Density (g/cc)	Balls
08295	Polyamide 6/6 (Nylon)	1/16"	1.19	100
19529	Polyethylene, Hollow	3/8"	0.53	100
16026	Polypropylene	1/4"	0.92	100
19827	Polystyrene, Smooth Finish	1/4"	1.05	100
17175	Polystyrene, Etched Surface	1/8"	1.05	100
16730	Polystyrene, Etched Surface	1/4"	1.05	100 or 500
18547	Polystyrene, Biological Applications	1/8"	1.05	100
18548	Polystyrene, Biological Applications	1/4"	1.05	100
19842	Polystyrene, Blue Dyed	1/4"	1.05	50
19841	Polystyrene, Carboxylate	1/8"	1.05	50
19840	Polystyrene, Carboxylate	1/4"	1.05	50
17649	Teflon	1/8"	2.14	100
17650	Teflon	1/4"	2.14	100

To Order

In The U.S. Call: 1-800-523-2575 • 215-343-6484

In The U.S. FAX: 1-800-343-3291 • 215-343-0214

In Germany Call: (49) 6221-765767

In Germany FAX: (49) 6221-764620

Order online anytime at www.polysciences.com.