

## BioMag<sup>®</sup> Cell Sorting References

- Adkins, B., A. Ghanei, K. Hamilton.** 1993. Developmental regulation of IL-4, IL-2, and IFN-gamma production by murine peripheral T lymphocytes. *J Immunol*, 151(12):6617-6626.
- Adkins, B., K. Hamilton.** 1992. Freshly isolated, murine neonatal T cells produce IL-4 in response to anti-CD3 stimulation. *J Immunol*, 149(11):3448-3455.
- Agematsu, K., T. Kobata, K. Sugita, G.J. Freeman, M.P. Beckmann, S.F. Schlossman, C. Morimoto.** 1994. Role of CD27 in T cell immune response. *J Immunol*, 153(4):1421-1429.
- Agematsu, K., T. Kobata, K. Sugita, T. Hirose, S.F. Schlossman, C. Morimoto.** 1995. Direct cellular communications between CD45R0 and CD45RA T cell subsets via CD27/CD70. *J Immunol*, 154(4):3627-3635.
- Alzona, M., H.M. Jack, R.I. Fisher, T.M. Ellis.** 1995. IL-12 activates IFN-g production through the preferential activation of CD30+ T cells. *J Immunol*, 154(1):9-16.
- Askenase, P.W., M. Szczepanik, M. Ptak, V. Paliwal, W. Ptak.** 1995. Gamma delta T cells in normal spleen assist immunized ab T cells in the adoptive cell transfer of contact sensitivity. *J Immunol*, 154(8):3644-3653.
- Bathon, J.M., J.C. Croghan, D.W. MacGlashan, D. Proud.** 1994. Bradykinin is a potent and relatively selective stimulus for cytosolic calcium elevation in human synovial cells. *J Immunol*, 153(6):2600-2608.
- Beno, D.W., A.G. Stover, H.L. Mathews.** 1995. Growth inhibition of *Candida albicans* hyphae by CD8+ lymphocytes. *J Immunol*, 154(10):5273-5281.
- Berman, J.S., D.M. Center.** 1987. Chemotactic activity of porcine insulin for human T lymphocytes in vitro. *J Immunol*, 138(7):2100-2103.
- Bieva, C.J., F.J. VanderBruggen, P.A. Stryckmans.** 1989. Malignant leukemic cell separation by iron colloidal immunomagnetic adsorption. *Experimental Hematology*, 17(8):914-920.
- Bonomo, A., P.J. Kehn, E. Payer, L. Rizzo, A.W. Cheever, E.M. Shevach.** 1995. Pathogenesis of post-thymectomy autoimmunity: role of syngeneic MLR-reactive T cells. *J Immunol*, 154(12):6602-6611.
- Brissette-Storkus, C. et al.** 1994. Characterization and function of the NKR-P1dim/T cell receptor-ab+ subset of rat T cells. *J Immunol*, 152(2):388-396.
- Chen, Y., V.K. Kuchroo, J. Inobe, D.A. Hafler, H.L. Weiner.** 1994. Regulatory T cell clones induced by oral tolerance: suppression of autoimmune encephalomyelitis. *Science*, 265(5176):1237-1240.
- Chervenak, R., D. Dempsey, R. Soloff, R.M. Wolcott, S.R. Jennings.** 1993. The expression of CD4 by T cell precursors resident in both the thymus and the bone marrow. *J Immunol*, 151(9):4486-4493.
- Chervonsky, A.V., T.V. Golovkina, S.R. Ross, C.A. Janeway.** 1995. Differences in the avidity of TRC interactions with a superantigenic ligand affect negative selection, but do not allow positive selection. *J Immunol*, 155(11):5115-5123.
- Chilton, P.M., R. Fernandez-Botran.** 1993. Production of soluble IL-4 receptors by murine spleen cells is regulated by T cell activation and IL-4. *J Immunol*, 151(11):5907-5917.
- Conlon, K.C., A.C. Ochoa, W.C. Kopp, J.R. Ortaldo, W.J. Urba, D.L. Longo, H.A. Young.** 1992. Enhanced lymphokine production and lymphokine receptor expression in multiple antibody-stimulated human CD4+ peripheral blood lymphocytes. *J Immunol*, 149(10):3278-3289.
- Cruikshank, W.W., J.S. Berman, A.C. Theodore, J. Bernardo, D.M. Center.** 1987. Lymphokine activation of T4+ T lymphocytes and monocytes. *J Immunol*, 138(11):3817-3823.
- Dieli, F., G.L. Asherson, G.C. Romano, G. Sireci, F. Gervasi, A. Salerno.** 1994. IL-4 is essential for the systemic transfer of delayed hypersensitivity by T cell lines. *J Immunol*, 152(6):2698-2704.
- Dong, Z., H.H. Wortis.** 1994. Function of bone marrow stromal cell lines derived from nude mice. *J Immunol*, 153(4):1441-1448.
- Elson, C.O., S.P. Holland, M.T. Dertzbaugh, C.F. Cuff, A.O. Anderson.** 1995. Morphologic and functional alterations of mucosal T cells by cholera toxin and its B subunit. *J Immunol*, 154(3):1032-1040.
- Ettinger, R. et al.** 1994. Functional distinctions between MRL-lpr and MRL-gld lymphocytes. *J Immunol*, 152(4):1557-1568.

- Ferbas, J.J., J.F. Toso, A.J. Logar, J.S. Navratil, C.R. Rinaldo.** 1994. CD4+ blood dendritic cells are potent producers of IFN-alpha in response to invitro HIV-1 infection. *J Immunol*, 152(9):4649-4662.
- Ford, J.R., M Terzaghi-Howe.** 1992. Characteristics of magnetically separated rat tracheal epithelial cell populations, *Am J Physiol*, 263:L568-L574.
- Fowler, D.H., K. Kurasawa, A. Husebekk, P.A. Cohen, R.E. Gress.** 1994. Cells of Th2 cytokine phenotype prevent LPS-induced lethality during murine graft-versus-host reaction. *J Immunol*, 152(3):1004-1013.
- Freedman, M.S., M. Blain, J.P. Antel.** 1991. Differential responses of CD4+CD45RA+ and CD4+CD29+ subsets to activated CD8+ cells: enhanced stimulation of the CD4+CD45RA+ subset by cells from patients with Multiple Sclerosis. *Cellular Immunology*, 133(2):306-316.
- Gilat, D., R. Hershkoviz, Y.A. Mekori, I. Vlodaysky, O. Lider.** 1994. Regulation of adhesion of CD4+ T lymphocytes to intact or heparinase-treated subendothelial extracellular matrix by diffusible or anchored RANTES and MIP-1b. *J Immunol*, 153(11):4899-4906.
- Gonzalez, M., R. Merino, A.L. Gonzalez, J. Merino.** 1995. The ability of B cells to participate in allogeneic cognate T-B cell interactions in vitro depends on the presence of CD4+ T cells during their development. *J Immunol*, 155(3):1091-1100.
- Gribben, J.G., et al.** 1992. Bone marrows of non-Hodgkin's lymphoma patients with a bcl-2 translocation can be purged of polymerase chain reaction-detectable lymphoma cells using monoclonal antibodies and immunomagnetic bead depletion. *Blood*, 80(4):1083-1089.
- Gryllis, C., M.A. Wainberg, M. Gornitsky, B. Brenner.** 1990. Diminution of inducible lymphokine-activated killer cell activity in individuals with AIDS-related disorders. *AIDS*, 4(12):1205-1212.
- Haregewoin, A., G. Soman, R.C. Horm, R.W. Finberg.** 1989. Human gamma delta+ T cells respond to mycobacterial heat shock protein. *Nature*, 340(6231):309-312.
- Hermentin, P., et al.** 1990. Hinge-thiol coupling of monoclonal antibody to silanized iron oxide particles and evaluation of magnetic cell depletion. *Bioconjugate Chemistry*, 1(6):411-418.
- Johnston, J.A., D.D. Taub, A.R. Lloyd, K. Conlon, J.J. Oppenheim, D.J. Kevlin.** 1994. Human T lymphocyte chemotaxis and adhesion induced by vasoactive intestinal peptide. *J Immunol*, 153(4):1762-1768.
- June, C.H., J.A. Ledbetter, M.M. Gillespie, T. Lindsten, C.B. Thompson.** 1987. T cell proliferation involving the CD28 pathway is associated with cyclosporine-resistant Interleukin-2 gene expression. *Molecular and Cell Biology*, 7(12):4472-4480.
- Kobata, T., K. Agematsu, J. Kameoka, S.F. Schlossman, C. Morimoto.** 1994. CD27 is a signal-transducing molecule involved in CD45RA+ naive T cell costimulation. *J Immunol*, 153(12):5422-5432.
- Kornfeld, H., W.W. Cruikshank, S.W. Pyle, J.S. Berman, D.M. Center.** 1988. Lymphocyte activation by HIV-1 envelope glycoprotein. *Nature*, 335(6189):445-448.
- Laberge, S., W.W. Cruikshank, D.J. Beer, D.M. Center.** 1996. Secretion of IL-16 (lymphocyte chemoattractant factor) from serotonin-stimulated CD8+ T cells in vitro. *J Immunol*, 156(1):310-315.
- Lagoo, A.S., et al.** 1994. IL-2, IL-4, & IFN-gamma gene expression versus secretion in superantigen-activated T cells. *J Immunol*, 152(4):1641-1652.
- Lederman, S., et al.** 1994. T-BAM/CD40-L on helper T lymphocytes augments lymphokine-induced B cell Ig isotype switch recombination and rescues B cells from programmed cell death. *J Immunol*, 152(5):2163-2171.
- Lee, W.T., E.S. Vitetta.** 1991. The differential expression of homing and adhesion molecules on virgin and memory T cells in the mouse. *Cellular Immunology*, 132(1):215-222.
- Levin, D., S. Constant, T. Pasqualini, R. Flavell, K. Bottomly.** 1993. Role of dendritic cells in the priming of CD4+ lymphocytes to peptide antigen in vivo. *J Immunol*, 151(12):6742-6750.
- Linette, G.P., R.J. Hartzman, J.A. Ledbetter, C.H. June.** 1988. HIV-1-infected T cells show a selective signalling defect after perturbation of CD3/antigen receptor. *Science*, 241(4865):573-576.
- Loeffler, C.M., et al.** 1991. Antitumor effects of interleukin 2 liposomes and anti-CD3-stimulated T cells against murine MCA-38 hepatic metastasis. *Cancer Res*, 51(8):2127-2132.
- Lopes, M.F., V.F. deVeiga, A.R. Santos, M.E. Fonseca, G.A. DosReis.** 1995. Activation-induced CD4+ T cell death by apoptosis in experimental Chagas' disease. *J Immunol*, 154(2):744-752.
- Mehrotra, P.T., A.J. Grant, J.P. Siegel.** 1995. Synergistic effects of IL-7 and IL-12 on human T cell activation. *J Immunol*, 154(10):5093-5102.
- Mokhtarian, F., Y. Shi, D. Shirazian, L. Morgante, A. Miller, D. Grob.** 1994. Defective production of anti-inflammatory cytokine, TGF-beta, by T cell lines of patients with active multiple sclerosis. *J Immunol*, 152(12):6003-6010.
- Moran, T.M., E. Toy, Y. Kuzu, H. Kuzu, H. Isobe, J.L. Schulman.** 1995. Staphylococcal enterotoxin B-activated T cells can be redirected to inhibit multicycle virus replication. *J Immunol*, 155(2):759-765.

- Nagata, M., P. Santamaria, T. Kawamura, T. Utsugi, J.W. Yoon.** 1994. Evidence for the role of CD8+ cytotoxic T cells in the destruction of pancreatic beta-cells in non-obese diabetic mice. *J Immunol*, 152(4):2042-2050.
- Nestle, F.O., X.G. Zheng, C.B. Thompson, L.A. Turka, B.J. Nickoloff.** 1993. Characterization of dermal dendritic cells obtained from normal human skin reveals phenotypic and functionally distinctive subsets. *J Immunol*, 151(11):6535-6545.
- Padmanabhan, R., C.D. Corsico, T.H. Howard, W. Holter, C.M. Fordis, M. Willingham, B.H. Howard.** 1989. Purification of transiently transfected cells by magnetic affinity cell sorting. *Analytical Biochemistry*, 170(2):341-348.
- Padmanabhan, R., C. Corsico, W. Holter, T.H. Howard, B.H. Howard.** 1989. Purification of transiently transfected cells by magnetic-affinity cell sorting. *Journal of Immunogenetics*, 116(2):91-102.
- Pilaro, A.M., D.D. Taub., K.L. McCormick, H.M. Williams, T.J. Sayers, W.E. Fogler, R.H. Wiltrot.** 1994. TNF-alpha is a principle cytokine involved in the recruitment of NK cells to liver parenchyma. *J Immunol*, 153(1):333-342.
- Pricop, L., et al.** 1993. Characterization of the Fc mu receptor on human natural killer cells. Interaction with its physiologic ligand, human normal IgM, specificity of binding, and functional effects. *J Immunol*, 151(6):3018-3029.
- Rabinowich, H., W.C. Lin, R.B. Herberman, T.L. Whiteside.** 1994. Signaling via CD7 molecules on human NK cells. Induction of tyrosine phosphorylation and beta/integrin-mediated adhesion to fibronectin. *J Immunol*, 153(8):3504-3513.
- Rabinowich, H., P. Sedlmayr, R.B. Herberman, T.L. Whiteside.** 1993. Response of human NK cells to IL-6 alterations of the cell surface phenotype, adhesion to fibronectin and laminin, and tumor necrosis factor-alpha/beta secretion. *J Immunol*, 150(11):4844-4855.
- Rabinowich, H., et al.** 1992. Role of cytokines in the adoptive immunotherapy of an experimental model of human head and neck cancer by human IL-2-activated natural killer cells. *J Immunol*, 149(1):340-349.
- Rabinowich, H., W.C. Lin, A. Amoscato, R.B. Herberman, T.L. Whiteside.** 1995. Expression of vitronectin receptor on human NK cells and its role in protein phosphorylation, cytokine production, and cell proliferation. *J Immunol*, 154(3):1124-1135.
- Rasmussen, R.A., S.L. Counts, J.F. Daley, S.F. Schlossman.** 1994. Isolation and characterization of CD6- T cells from peripheral blood. *J Immunol*, 152(2):527-536.
- Rosenman, S.J., A.A. Ganji, T.F. Tedder, W.M. Gallatin.** 1993. Syn-capping of human T lymphocyte adhesion/activation molecules and their redistribution during interaction with endothelial cells. *J Leuk Biol*, 53(1):1-10.
- Saha, K., P.K. Wong.** 1992. Ts1, a temperature-sensitive mutant of Moloney murine leukemia virus TB, can infect both CD4+ and CD8+ T cells, but requires CD4+ T cells in order to cause paralysis and immunodeficiency. *J Virol*, 66(5):2639-2646.
- Saitoh, S.I., A. Kosugi, S. Noda, N. Yamamoto, M. Ogata, Y. Minami, K. Miyake, T. Hamaoka.** 1995. Modulation of TCR-mediated signaling pathway by thymic shared antigen-1 (TSA-1)/stem cell antigen-2 (Sca-2). *J Immunol*, 155(12):5574-5581.
- Santamaria, P., T. Utsugi, B. Park, N. Averill, S. Kawazu, J.W. Yoon.** 1995. Beta-cell-cytotoxic CD8+ T cells from non-obese diabetic mice use highly homologous T cells receptor alpha-chain CDR3 sequences. *J Immunol*, 154(5):2494-2503.
- Sen-Majumdar, A., C. Guidos, T. Kina, M. Lieberman, I.L. Weissman.** 1994. Characterization of preneoplastic thymocytes and of their neoplastic progression in irradiated C57BL/Ka mice. *J Immunol*, 153(4):1581-1592.
- Selvey, L.A., H.C. Morse, C.H. June, R.J. Hodes.** 1995. Analysis of antigen receptor signalling in B cells from mice with retrovirus-induced acquired immunodeficiency syndrome. *J Immunol*, 154(1):171-179.
- Shen, F., X.L. Xu, L.H. Graf, S.F. Chong.** 1995. CD45-cross-linking stimulates IFN-gamma production in NK cells. *J Immunol*, 154(1):644-652.
- Smyth, M.J., S.L. Strobl, H.A. Young, J.R. Ortaldo, A.C. Ochoa.** 1991. Regulation of lymphokine-activated killer activity and pore forming protein gene expression in human peripheral blood CD8+ T lymphocytes. Inhibition by transforming growth factor-beta. *J Immunol*, 146(10):3289-3297.
- Spangrude, G.J., R. Scollay.** 1990. A simplified method for enrichment of mouse hematopoietic stem cells. *Experimental Hematology*, 18(8):920-926.
- Spangrude, G.J., S. Heimfeld, I.L. Weissman.** 1988. Purification and characterization of mouse hematopoietic stem cells. *Science*, 241(4861):58-62.
- Stout, R.D., J. Suttles, J. Xu, I.S. Grewal, R.A. Flavell.** 1996. Impaired T cell-mediated macrophage activation in CD40 ligand-deficient mice. *J Immunol*, 156(1):8-11.
- Sugie, K., Y. Minami, T. Kawakami, A. Uchida.** 1995. Stimulation of NK-like YT cells via leukocyte function-associated antigen (LFA)-1. Possible involvement of LFA-1-associated tyrosine kinase in signal transduction after recognition of NK target cells. *J Immunol*, 154(4):1691-1698.
- Tacke, M., G.J. Clark., M.J. Dallman, T. Hunig.** 1995. Cellular distribution and costimulatory function of rat CD28. Regulated expression during thymocyte maturation and induction of cyclosporin A sensitivity of costimulated T cell responses by phorbol ester. *J Immunol*, 154(10):5121-5127.
- Taylor, R.P., F. Pocanic, C. Reist, E.L. Wright.** 1991. Complement-opsonized IgG antibody/dsDNA immune complexes bind to CR1 clusters on isolated human erythrocytes. *Clinical Immunology and Immunopathology*, 61(2 Pt. 1):143-160.

Should any of our materials fail to perform to our specifications, we will be pleased to provide replacements or return the purchase price. We solicit your inquiries concerning all needs for life sciences work. The information given in this bulletin is to the best of our knowledge accurate, but no warranty is expressed or implied. It is the user's responsibility to determine the suitability for their own use of the products described herein, and since conditions of use are beyond our control, we disclaim all liability with respect to the use of any material supplied by us. Nothing contained herein shall be construed as a recommendation to use any product or to practice any process in violation of any law or any government regulation.

- Trickett, A.E., D.J. Ford, P.R. Lam-Po-Tang, M.R. Vowels.** 1990. Comparison of magnetic particles for immunomagnetic bone marrow purging using an acute lymphoblastic leukemia model. *Transplantation Proceedings*, 22(5):2177-2178.
- Turka, L.A., D.G. Schatz, M.A. Oettinger, J.J. Chun, C. Gorka, K. Lee, W.T. McCormack, C.B. Thompson.** 1991. Thymocyte expression of RAG-1 and RAG-2: termination by T cell receptor cross-linking. *Science*, 253(5021):778-781.
- Vaccaro, D.E.** 1990. Applications of magnetic separation: cell sorting. *Am BiotechnolLab*, 8(5):30, 32-35.
- Vujanovic, N.L., S. Yasumura, H. Hirabayashi, W.C. Lin, S. Watkins, R.B. Herberman, T.L. Whiteside.** 1995. Antitumor activities of subsets of human IL-2-activated natural killer cells in solid tissues. *J Immunol*, 154(1):281-289.
- Wirth, S., L.G. Guidotti, K. Ando, H.J. Schlicht, F.V. Chisari.** 1995. Breaking tolerance leads to autoantibody production, but not autoimmune disease in Hepatitis B virus envelope transgenic mice. *J Immunol*, 154(5):2504-2515.
- Witkowski, J.M., S.P. Li, G. Gorgas, R.A. Miller.** 1994. Extrusion of the P glycoprotein substrate rhodamine-123 distinguishes CD4 memory T cell subsets that differ in IL-2-driven IL-4 production. *J Immunol*, 153(2):658-665.
- Yabkowitz, R., V.M. Dixit, N. Guo, D.D. Roberts, Y. Shimizu.** 1993. Activated T cell adhesion to thrombospondin is mediated by the alpha 4 beta 1 (VLA-1) and alpha 5 beta 1 (VLA-5) integrins. *J Immunol*, 151(1):149-158.
- Yuen, M.H., M.P. Protti, B. Diethelm-Okita, L. Moiola, J.F. Howard, B.M. Conti-Fine.** 1995. Immunoregulatory CD8+ cells recognize antigen-activated CD4+ cells in myasthenia gravis patients and in healthy controls. *J Immunol*, 154(3):1508-1520.