

Hollow Glass Beads


Catalog Number: 19823

Description


Our hollow glass beads offer the strength and stability of glass at 1.1 g/cc nominal density.

Characteristics

Composition:	Proprietary glass oxide, non-porous
Mean Diameter:	8 μ m (2-20 μ m range)
Shape:	Spherical
Appearance:	White powder
Specific Gravity:	1.1 g/cc
Crush Strength:	Survives injection molding, >10,000 psi
Chemical Resistance:	Low alkali leach, insoluble in water



Hollow Glass Spheres 2-20 μ m at 100x magnification



Hollow Glass Spheres 2-20 μ m at 100x magnification

Storage and Stability

Store dry product at room temperature.

This product is for research use only and is not intended for use in humans or for *in vitro* diagnostic use.

Ordering Information

Cat. #	Description	Size
19823-5	Hollow Glass Beads	5g

To Order

In The U.S. Call: 1-800-523-2575 • 215-343-6484

In The U.S. FAX: 1-800-343-3291 • 215-343-0214

In Germany Call: (49) 6221-765767

In Germany FAX: (49) 6221-764620

Order online anytime at www.polysciences.com.

Should any of our materials fail to perform to our specifications, we will be pleased to provide replacements or return the purchase price. We solicit your inquiries concerning all needs for life sciences work. The information given in this bulletin is to the best of our knowledge accurate, but no warranty is expressed or implied. It is the user's responsibility to determine the suitability for their own use of the products described herein, and since conditions of use are beyond our control, we disclaim all liability with respect to the use of any material supplied by us. Nothing contained herein shall be construed as a recommendation to use any product or to practice any process in violation of any law or any government regulation.