



U.S. Corporate Headquarters
400 Valley Road
Warrington, PA 18976
1(800) 523-2575 / (215) 343-6484
Fax: 1(800)343-3291
Email: info@polysciences.com

Polysciences Europe GmbH
Handelsstrasse 3 D-69214
Eppelheim, Germany
+(49) 6221-765767
Fax: +(49) 6221-764620
Email: info@polysciences.de

TECHNICAL DATA SHEET 891

Page 1 of 1

High Temperature Resistant 1 Part Epoxy EW8002 N

EW 8002 N is a high temperature resistant epoxy. The product is well suited for many bonding and sealing applications requiring high modulus and high temperature resistance.

PRODUCT FEATURES

- High Thermal Resistance
- Rigid and High Strength
- Durable
- High Opacity

TYPICAL PROPERTIES

COMPOSITION	Epoxy
APPEARANCE	Black Liquid
VISCOSITY @25° C	30,000 cps
SPECIFIC GRAVITY	1.70

OPTIMUM CURE CYCLE

90 minutes @ 100° C + 60 minutes @ 165° C

CURED PROPERTIES

TG	150° C
CTE	25 ppm

USAGE

EW8002 N is stored frozen to extend its shelf life, the product must be thawed prior to use. Apply the thawed product to the bonding surface and cure at the recommended curing cycle.

SHELF LIFE

6 months frozen

TO ORDER

In The U.S. Call:	1(800) 523-2575 • (215) 343-6484
In The U.S. Fax:	1(800) 343-3291 • (215) 343-0214
In Germany Call:	+(49) 6221-765767
In Germany Fax:	+(49) 6221-764620