



U.S. Corporate Headquarters  
400 Valley Rd.  
Warrington, PA 18976  
1(800) 523-2575 / (215) 343-6484  
1(800)343-3291 fax  
info@polysciences.com

Polysciences Europe GmbH  
Handelsstrasse 3  
D-69214 Eppelheim, Germany  
+(49) 6221-765767  
+(49) 6221-764620 fax  
info@polysciences.de

Polysciences Asia-Pacific, Inc.  
2F-1, 207 Tunhwa N. Rd.  
Taipei, Taiwan 10595  
(886) 2 8712 0600  
(886) 2 8712 2677 fax  
info@polysciences.tw

*For Immediate Release*

## **PRESS RELEASE**

1/01/2014

Contact Information:  
**Jeff DePinto**  
Vice President Catalog  
& Intl. Business Development  
[jdepinto@polysciences.com](mailto:jdepinto@polysciences.com)

## **Polysciences, Inc. Releases 2014 – 2015 Product Catalog**

***Featuring hundreds of new products!***

**WARRINGTON, PA, Jan, 2014** – Polysciences, Inc., a leading manufacturer of laboratory supplies and specialty monomers and polymers for biotechnology, medical research and analytical chemistry industries, has published its new 456 page full-color 2014-2015 Product Catalog. This new catalog represents four product families and introduces hundreds of new product offerings for Life Sciences, Monomers & Polymers, Microspheres & Particles and High Performance Adhesives markets. The 2014-2015 Catalog lists thousands of off-the-shelf products that are used by researchers, lab technicians, chemists, scientists and engineers in the fields of histology, microscopy, life sciences, pharmaceuticals and polymer chemistry. These products include fine and specialty chemicals, fluorescent and certified dyes, monomers and polymers, magnetic and non-magnetic particles, microspheres & particle kits, and high performance adhesives & coatings.

### **What's New in the 2014-2015 Catalog!**

**Life Sciences** – Chemzymes Ultra Pure<sup>®</sup> Enzymes, Cell Culture Bags and Bovine Serum Albumin (BSA)

**Monomers & Polymers** – Expanded line of Hydrophilic Monomers, Ion Exchange Resins and Enzyme Carrier Resins for Enzyme Immobilization

**Microspheres & Particles** – ProMag<sup>®</sup> HP High Performance Magnetic Microspheres, Fluoresbrite Europium Chelate Carboxylate Microspheres, Small Bead Calibration Kits

**High Performance Adhesives & Encapsulants** – PC 220 Epoxy Encapsulant Kit and PC 7050 UV Curable Epoxy

To request a free copy, please visit [www.polysciences.com/newcat](http://www.polysciences.com/newcat) or call **1-800-523-2575**

### ***About Polysciences, Inc.***

*Polysciences Inc. was founded in 1961 and specializes in the development and manufacturing of high purity chemicals and reagents for laboratory research, histology, healthcare diagnostics, and medical devices. Their work in the biomedical and pharmaceutical industries has helped bring about countless innovations in healthcare, and their offerings in high performance adhesives and personal care products offer even more innovative opportunities for the company.*

*For more information about Polysciences, Inc., visit [www.polysciences.com](http://www.polysciences.com)*

###