

# Polystyrene Beads, Large


Catalog Numbers: 19825, 19826, 21392

## DESCRIPTION

Our crosslinked Polystyrene Beads are larger than Polybeads® (50nm – 90µm), yet smaller than Polyballs (1/16" – 1/4"). As non-functionalized polystyrene spheres, they may be coated with proteins using adsorption protocols for applications that call for large beads. Diameter ranges are approximate and represent sieved fractions.

## CHARACTERISTICS

Size Range 106 – 600µm



Cat. #21392: Polystyrene Beads, Large (500-600µm)

## STORAGE

Store at room temperature.

This product is for research use only and is not intended for use in humans or for *in vitro* diagnostic use.

## ORDERING INFORMATION

Cat. #	Description	Size
19825	Polystyrene Beads, Large (200-300µm-nominal)	1g
19826	Polystyrene Beads, Large (355-425µm-nominal)	1g
21392	Polystyrene Beads, Large (500-600µm-nominal)	1g

Visit [Polysciences.com](https://www.polysciences.com) anytime to place an order.

Should any of our materials fail to perform to our specifications, we will be pleased to provide replacements or return the purchase price. We solicit your inquiries concerning all needs for life sciences work. The information given in this bulletin is to the best of our knowledge accurate, but no warranty is expressed or implied. It is the user's responsibility to determine the suitability for their own use of the products described herein, and since conditions of use are beyond our control, we disclaim all liability with respect to the use of any material supplied by us. Nothing contained herein shall be construed as a recommendation to use any product or to practice any process in violation of any law or any government regulation.