

Total Cutting Solution Kit



Cut Difficult Tissue with Ease!

A revolutionary product for every work station in your lab. Our **Total Cutting Solution Kit** provides you with 5 tissue cutting solutions to address decalcification issues and difficult tissue blocks & specimens. Don't let poor processing and paraffin issues stop you in the middle of your cutting track record. Reduce bottlenecks in your workflow by using the 10 tissue block wells provided with the kit right at your cutting station!

The Total Cutting Solution Kit is the answer for all work stations!



Histoheme, Soft Block & Soft Nail

- Soften paraffin blocks that are dry, brittle and hard
- Hematological samples cut in long clean ribbons
- Reduces faded appearance of stains
- Break down structural protein of hair, nail, hoof, feathers and skin
- Great for Mohs Labs, Veterinary Histology & Podiatry Nail Histology

Super Decalcifier I: Delicate

- Recommended for use with delicate specimens that require fast diagnostic work such as IHC, cytochemistry and special stain procedures.

Super Decalcifier II: Heavy Duty

- Designed for safe, longer term decalcification of hard compact bone, femoral heads, above and below the knee amputations.

Cat. #: 25388

Size: 1 kit

We offer a complete line of BioSciences products, to learn more visit polysciences.com