

Flotation Work Stations

Make your tissue preparation easier, safer and more efficient. Removable glass dish waterbath with chemical resistant plastic housing, optimum-specimen viewing with a high contrast background illuminated by an LED light array. The dependability and safety of up-to-date microprocessor controlled electronics, combined with histo/orientor for eliminating wrinkles in sections and a heated slide dryer in one compact unit.

Flotation Work Station with 8" x 8" x 2 1/4" (Deep Dish), HistoOrientor and Dryer

Cat. # 25389-1



Flotation Work Station, microprocessor controlled with removable glass bowl w/built-in LED light array, 2" x 5" slide dryer and 2" x 2" HISTO/Orientator™, TUV/CE/FCC.

Specifications: 100-120VAC ± 10%, 50/60Hz, 3.7A

Dimensions (L x W x H): 46.99 cm x 30.73 cm x 14.22 cm

Weight: 3.99 kg (8.8 lbs)

continued on back

Flotation Work Stations

Make your tissue preparation easier, safer and more efficient. Removable glass dish waterbath with chemical resistant plastic housing, optimum-specimen viewing with a high contrast background illuminated by an LED light array. The dependability and safety of up-to-date microprocessor controlled electronics, combined with histo/orientor for eliminating wrinkles in sections and a heated slide dryer in one compact unit.

Flotation Work Station with 8" x 11" x 2" Dish HistoOrientor and Dryer

Cat. # 26011-1

Specifications: 100-120VAC \pm 10%, 50/60Hz, 5.4A

Dimensions (L x W x H): 47.62 cm x 38.1 cm x 13.97 cm

Weight: 4.63 kg (10.2 lbs)



Flotation Work Station with 8" x 11" x 2" Dish

Cat. # 26010-1

Specifications: 100-120 VAC \pm 10%, 50/60Hz, 5A

Dimensions (L x W x H): 47.62 cm x 38.1 cm x 13.97 cm

Weight: 4.15 kg (9.15 lbs)



Flotation Work Station with 8" x 8" x 2 1/4" (Deep Dish)

Cat. # 26012-1

Specifications: 100-120 VAC \pm 10%, 50/60Hz, 3.4A

Dimensions (L x W x H): 46.99 cm x 30.73 cm x 14.22 cm

Weight: 3.49 kg (7.7 lbs)



Flotation Work Station, microprocessor controlled with removable glass bowl w/built-in LED light array, 2" x 5" slide dryer and 2" x 2" HISTO/Orientator™, TUV/CE/FCC.