


Paraffin

GemCut™ Paraffins

Cat. #24361-24366

GemCut™ Paraffins, with a melting point range of 56-58°C, is a new high polymer paraffin formulation for routine embedding and cutting. This paraffin is specially formulated to exhibit superior cutting performance in high humidity environments, with low compression. The paraffin is available in multiple colors to assist the histologic technician in visualizing small biopsies. The coloring agents and polymer content do not interfere with staining protocols.



Available in six color formulations for visualizing small biopsies.

[Technical Data Sheet #691](#)

Cat. #	Description	Size
24361-1	GemCut™ Pink Sapphire	1 case (4x2kg)
24363-1	GemCut™ Citrine	1 case (4x2kg)
24364-1	GemCut™ Emerald	1 case (4x2kg)

PolyFin® Paraffin

Cat. #19562

A mixture of fine paraffin waxes and copolymer alloys in convenient pellets. PolyFin® produces ultra-thin compression-free ribbons. It is used to cut 2 micron sections routinely without icing. It functions flawlessly in pressure/vacuum fluid flow processors. PolyFin® provides maximum support for both hard and soft tissue samples while maintaining exceptional clarity. The low 55°C melting point eliminates excessive heat distortion during tissue processing. [Technical Data Sheet #429](#)



Cat. #	Description	Size
19562	PolyFin® Paraffin	1 case (8x1kg) / 15kg

Peel-A-Way® Micro-Cut Paraffins

Cat. #24201, #24198, #24202

Micro-Cut paraffin, available at 3 melting temperatures is a low polymer formulation in our Peel-A-Way® paraffin line. This low polymer paraffin gives improved viscosity, infiltration and embedding. The translucent appearance of the paraffin allows excellent tissue visibility and easier paraffin removal with all types of clearing agents. Available in a larger, more convenient 10kg (22lb) case size containing four 2.5kg stand alone bags. One case will fill most tissue processors. Available in 3 melting points. [Technical Data Sheet #525A](#)



Cat. #	Description	Size
24201-1	Micro-Cut Paraffin 52°-54°	1 case (4x2.5kg)
24362-1	Micro-Cut Paraffin 56°-58°	1 case (4x2.5kg)
24363-1	Micro-Cut Paraffin 62°-64°	1 case (4x2.5kg)

Peel-A-Way® Micro-Cut Preferred (DMF) Paraffin

Cat. #24198D

Micro-Cut Preferred paraffin, is a low polymer formulation with DMF (dimethylformamide) added to increase infiltration of tissue.



Cat. #	Description	Size
24198D-1	Micro-Cut Preferred 56°-58°	1 can (4x2.5kg)

Embedding Mold Release Spray

Cat. #03928

A particularly effective anti-stick spray that dries on contact and leaves no oily film, allowing quick release of specimens from molds. It is unaffected by water, oil or solvents.



Cat. #	Description	Size
03928-1	Mold Release Spray	1 can