



Phen-DC3

G-Quadruplex Ligand

**Polysciences, Inc. is a new commercial supplier of Phen-DC3.
Please ask us about synthesizing additional G4 ligands.**

Polysciences, Inc. provides custom manufacturing and development of fine chemicals and intermediates for the drug discovery, biotechnology, medical device, cosmetics, electronics and analytical chemistry industries. In addition, Polysciences provides custom synthesis solutions for research and development compounds on the milligram to ton scale. Our FDA audited production facilities are capable of both Good Laboratory Practice (GLP) and c-GMP level manufacturing.

polysciences.com