

StarLight™ Calibration Slides



Our StarLight™ Calibration Slides feature **vibrant ~6µm fluorescent microspheres** dyed with a single fluorophore for basic imaging checks and calibrations. Our **four standard versions** are appropriate for use with common microscope filter sets: **Glacial Blue** (360, 450), **Dragon Green** (480, 520), **Envy Green** (525, 565), **Flash Red** (660, 690), and are available individually or as the full 4-slide collection.








Figure 1: Images of StarLight™ Calibration Slides. The left image is Glacial Blue. The right images are (from top to bottom) Dragon Green, Envy Green, and Flash Red.

SPECTRA INFORMATION



For additional information regarding our accessory reagents, visit our website at www.polysciences.com.

StarLight™ Calibration Slides

Cat. #	Description
25441	StarLight™ Calibration Slide • Glacial Blue
25442	StarLight™ Calibration Slide • Dragon Green
25443	StarLight™ Calibration Slide • Envy Green
25444	StarLight™ Calibration Slide • Flash Red
25445	StarLight™ Collection • 4-Slide Pack

• StarLight™ is a trademark of Bangs Laboratories, Inc.

Order online anytime at www.polysciences.com