

## TECHNICAL DATA SHEET 989

Page 1 of 1

# 3-Phenoxy 2 hydroxy propyl methacrylate (PHPM)

CAS #: 16926-87-7

### APPLICATIONS

Hydrophilic monomer useful in medical device and ophthalmic applications.

### CHEMICAL FORMULA

$C_{13}H_{16}O_4$


### MOLECULAR WEIGHT

236.26

### SYNONYMS

2-Hydroxy-3-Phenoxypropyl Methacrylate

2-Hydroxy-3-Phenoxypropylester



### ANALYSIS

%PHPM  
APHA color report  
% H<sub>2</sub>O  
pH  
MEHQ  
BHT  
Polymer  
Visual

### SPECIFICATIONS

95.0% min  
Less than 100 typical  
0.7% max  
6.0 - 8.0  
2500 ppm max  
1000 ppm max  
Presence  
Light yellow to slightly amber liquid

### ORDERING INFORMATION

Cat. #	Description	Size
25506	3-Phenoxy 2 hydroxy propyl methacrylate	100 g 500 g

### TO ORDER

In The U.S. Call:	1(800) 523-2575 • (215) 343-6484
In The U.S. Fax:	1(800) 343-3291 • (215) 343-0214
In Germany Call:	+(49) 6221-765767
In Germany Fax:	+(49) 6221-764620
In Asia Call:	(886) 2 8712 0600
In Asia Fax:	(886) 2 8712 2677

Order online anytime at [www.polysciences.com](http://www.polysciences.com)