



Leverage our facilities and expertise!

PROCESSING

- Solids & liquid filling
- Room temp & hot-pour, emulsion and gel materials
- Flammable & hazardous
- WFI grade and DI water systems
- Custom blending – solids, liquids, high viscosity
- Polymer processing (cryogenic milling)

PACKAGING

- Shrink tunnels & induction sealers, including tamper evident packaging
- Bulk drums, bottles, tubes, & packs
- In-line labeling
- Kit assembly

Through our Contract Manufacturing services, Polysciences offers tailored solutions to support your product manufacturing pipeline and accelerate time to market. Whether you're embarking on a new product launch, looking to us as a new source for current products or expanding your catalog through private labeling, we have the expertise, capabilities and operations to get you there. Our comprehensive and flexible services allow you to select exactly the service you need at any point in your workflow, and our state-of-the-art manufacturing and processing equipment can accommodate your project from pilot scale to bulk production.

As dedicated partners, our experienced project management team will be with you throughout your entire process. We're happy to see your project through formulation and packaging challenges, filling with options including world class ISO 5 or 8 cleanrooms, technical concerns, product design and more, giving our full attention to the project at hand, both quickly and accurately. Thanks to decades of experience, numerous clients big and small have successfully started, completed and continued their projects with us here at Polysciences.

When it comes to working with us, know that your project will be met with knowledge, ingenuity and care by our entire Contract Manufacturing team. With modern, flexible project spaces and production suites, cGMP and ISO 13485 certification, we're devoted to helping you develop superior product, all while meeting deadlines and bringing out the best in your concept.

From raw materials to the final product, we're prepared to be with you every step of the way.

