


Differential Quik Stain Kit (Modified Giemsa)



Cat.#: 24606



High quality, rapid turnaround time staining kit that serves many purposes in the cost conscious laboratory environment. The **Differential Quik Stain Kit** may be used for fine needle aspirations, cytology, blood smears, bone marrow biopsies as well as bone marrow smears, touch preps, various microorganisms, as well as detecting the *H. pylori* microorganism. The kit contains a fixative for air-dried cell suspensions and/or touch preparation slides.

Advantages:

- Interpret results in 15 seconds
- Visualize nuclear and cytoplasmic staining easily
- Differential stain designed to show virtually all dermapathology cells
- Rapid stain for immediate interpretation of fine needle aspiration biopsies
- Optimal as a hematology stain
- Excellent cytologic detail
- Consistently reliable results
- Detect blood-borne malarial parasites
- Available in 3 convenient kit sizes 250ml, 500ml and 1 gallon