


Tissue-Guard™ Coverslipping Film

*Save time in your day
by using Tissue-Guard;
the smart choice for your
coverslipping needs!*



Tissue-Guard Coverslipping Film provides efficiency in the high pace laboratory setting. The special designed plastic film provides exceptional clarity while protecting your specimen for long term storage. It is a safer alternative to glass coverslips. No breaking of coverslips or constant refilling of coverslip needed. Each roll of coverslipping film is 70 meters and can coverslip up to 1,200 slides. Tissue-Guard is compatible with all automated coverslippers requiring plastic film tape. For *In Vitro* Diagnostic use.

Catalog #25737-1

Size: 1 roll

Catalog #25737-5

Size: 1 case (5 rolls)