




BREAK THE MOLD

with Peel-A-Way® Disposable Embedding Molds




EMBEDDING MOLDS


Polysciences' Peel-A-Way Embedding Molds are useful for plastic, frozen, and epoxy embedding. The unique polyethylene construction is designed to "Peel-A-Way" from the solidified block, eliminating the need for release agents and block trimming.




TRUNCATED T8




TRUNCATED T12



SQUARE S22



RECTANGULAR R30



RECTANGULAR R40

ADVANTAGES

- Easy to use and saves time
- Unique design will "Peel-A-Way" from a solidified embedded block
- Versatile - use with PolyFreeze, JB-4®, GMA, MMA and epoxy.

Peel-A-Way Embedding Molds

Cat. #	Description
18985	Truncated - T8, 22mm x 22mm truncated to 8mm x 8mm, 20 mm deep
19896	Truncated - T12, 22mm x 22mm truncated to 12mm x 12mm, 20mm deep
18646A	Square - S22, 22mm x 22mm square x 20mm deep
18646B	Rectangular - R30, 22mm wide x 30mm long x 20mm deep
18646C	Rectangular - R40, 22mm wide x 40mm long x 20mm deep
18646D	Sampler Pack - Contains T8 (8), T12 (8), S22 (16), R30 (8), R40 (8)

Order today at [Polysciences.com](https://www.polysciences.com)