

BioMag® Multi-6 Microcentrifuge Tube Separator


Catalog Number: 8MB4111S

DESCRIPTION

The BioMag® Multi-6 Microcentrifuge Tube Separator is a durable, plastic magnetic separation unit designed to accommodate six 1.5ml microcentrifuge tubes. The magnets and tube holders are positioned to efficiently magnetically separate microliter volumes common in molecular biology applications.

PROCEDURE

Magnetic separations can be performed on six 1.5ml microcentrifuge tubes simultaneously. Once the tubes are placed in the tube holders, the magnetic particles are pulled to the side of the unit to complete a separation in less than one minute. Upon completion of the separation, the supernatant may be aspirated, leaving the magnetic particles as a pellet on the side of the tube.



STORAGE AND SAFETY

Storage The BioMag® Multi-6 Microcentrifuge Tube Separator should be wiped clean with mild detergent. Do not immerse. Harsh chemicals and organic solvents may damage the surface of the unit.

Safety Strong magnetic field. This separation unit should be used and stored a safe distance from magnetic media, such as computers, videotapes, credit cards, etc.

This product is for research use only and is not intended for use in humans or for *in vitro* diagnostic use.

ORDERING INFORMATION

Cat. #	Description	Size
8MB4111S-1	BioMag® Multi-6 Microcentrifuge Tube Separator	1 each

TO ORDER

In The U.S. Call: 1(800) 523-2575 • (215) 343-6484

In The U.S. Fax: 1(800) 343-3291 • (215) 343-0214

In Germany Call: +(49) 06201-845200

In Germany Fax: +(49) 06201-8452020

In Asia Call: (886) 2 8712 0600

In Asia Fax: (886) 2 8712 2677

Order online anytime at www.polysciences.com